


GLUCOSE MONITORS FOR DAILY DIABETES TESTING

Advanced technology glucose monitors kits or meters.
 Plus model is available as unit only or complete kit.
 EZ model only in complete kit. Kit includes 15 strips, code chip, control solution, 15 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manuals, logbook.

Main features:

- Accurate result in just 10 seconds.
- Only 1 µl sample size.
- Automatic calibration with code chip.
- Plasma calibration results.
- Advanced data management with 7, 14, 30 days average (only plus model).
- Affordable for frequent monitoring.
- Flexibility to test on the forearm, palm or fingertip.

- BIOSENSOR TECHNOLOGY
 - 10 SECOND TEST TIME
 - ONLY 1 µL BLOOD SAMPLE



23891

23900

		TECHNICAL SPECIFICATIONS	
Measurement Range		20 to 600 mg/dL or 1.1 to 33.3 mmol/L	
Sample Sites		Finger tips, Forearm or palm	
Sample Size		1 µL	
Test Time		Approximately 10 seconds	
Power Source		One 2032 or equivalent 3.0V coin cell battery	
Battery Life		12 months or approximately 1,000 tests	
Memory	Plus	Up to 300 records with time and date	
	EZ	Last 30 results, no averaging, no date/time	

				mg/dL	mmol/L	ON CALL® GLUCOSE MONITOR
Italian	English	Spanish	English			
Greek	French	Portuguese	French			
23900	23899	-	-			Glucose monitor EZ-meter only
23901*	-	-	-			Glucose monitor EZ-kit
23890	23891	23892	23881			Glucose monitor Plus-meter only
23895	23896	23897	23886			Glucose monitor Plus-kit

* Sold in Italian/Greek box with 6 languages manual (Italian, Greek, English, French, Spanish, Portuguese)

GIMA code	Accessories and consumable for GLUCOSE MONITOR
23910	Glucose Strips - box of 25 strips
23912	Glucose Strips - box of 50 strips
23913	Glucose Strips - box of 100 strips
23914	Control Solution
23915	Auto lancet device "On Call"
23916	Lancets 28G for 23915 - box of 100 pcs.



23896

SAFETY AUTOMATIC LANCETS

ONE-STEP SAFETY AUTOMATIC LANCETS

Multi purpose capillary blood sampling devices for all testing requirements. The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable. High-speed delivery and penetration method minimizes patient pain and operator error. Latex-free, hypo-allergenic, disposable.



24515

ONE-STEP SAFETY AUTOMATIC LANCETS							
GIMA code	Type	Depth of penetration	Gauge	Blood Flow	Finger Stick	Glucose	Hematocrit
24513	needle	1.0 mm	G26	5-10 µl	•	•	•
24515	needle	2.2 mm	G22	20-40 µl	•	•	•

LANCETS AND LANCET DEVICE

- 23915 AUTO LANCET DEVICE "ON CALL"
 Reusable auto lancet devices, in plastic for years use. Adjustable tip, allows precise configuration of lancet penetration.



23915

- 23940 LANCETS - 28G - box of 25 pcs - expires 12/2011
 - 23916 LANCETS - 28G - box of 100 pcs
- Lancets can be used with code 23915 and most of lancets devices in the market (see compatibility list in our website).



DEVICES AGAINST INSECTS, MEDUSAL, SMOKYBITES

• 25771 ZANZA-CLICK AFTERBITE

A few clicks to soothe the irritation caused by insects, medusae and nettle.



• 25772 BABYFRIEND

Ultrasonic repelling device against mosquitos.



• 25773 ECOSAVE

Electronic stimulator against snake bites or poisonous stings of insects and fishes. Few electric shocks are enough to reduce toxic reaction.



"MULTICARE IN" - SELF CONTROL OF RISK FACTOR


• 23965 MULTICARE "IN" KIT

A new approach for the periodical self-monitor of risk factors - cholesterol, triglycerides and glucose. The only device on sale with reflect-metrical and electro-chemical bio sensor suction. "Multicare in" is supplied in carton box with a nylon bag 10 glucose strips, lancet device, 10 lancets, rapid guide and user manual (English and Italian).

The first and only multi-parameter device with double technology


BIOSENSOR SUCTION STRIP - glucose

- Electro-chemical technology for a more practical use.
- Requires very little blood samples (1 µl).
- Short time for the results (5 sec.)



REFLECTOMETRIC SUCTION STRIP - cholesterol and tryglicerides

- multiCare in is a very singular technology.
- It allows a more functional self-monitoring of triglycerides and cholesterol.




3 IN 1: CHOLESTEROL, TRIGLYCERIDES, GLUCOSE BIOSENSOR



GIMA code	SPARE STRIPS FOR MULTICARE IN
23969	Glucose strips biosensor - tube of 25 strips + 1 chip
23970	Glucose strips biosensor - tube of 50 strips + 1 chip
23971	Cholesterol strips - tube of 5 strips + 1 chip
23972	Cholesterol strips - tube of 25 strips + 1 chip
23973	Triglycerides strips - tube of 5 strips + 1 chip
23974	Triglycerides strips - tube of 25 strips + 1 chip
LANCETS FOR MULTICARE AND MULTICARE IN	
23939	Multicare auto lancet device
23940	Lancets 28G - bag of 25 pcs - expiry 12/2011
23916	Lancets 28G - box of 100

TECHNICAL SPECIFICATIONS			
	Measuring range	Measuring time	Blood sample
Glucose electro-chemical:	10-600mg/dl	5"	1 µl
Cholesterol	150-400 mg/dl	30"	10 µl
Triglycerides	50-500 mg dl	30"	10 µl
Memory:	500 test results with date and time		
Batteries:	2 lithium batteries Cr 2032 3 Volts (32212)		
Size:	97x49x20.5 mm		
Weight:	65 g		

CHOLESTEROL, TRIGLYCERIDES, GLUCOSE POCKET MONITORS

- 23948 MULTICARE KIT - English/Spanish, manual
- 23949 MULTICARE KIT - Italian manual
- 23950 MULTICARE KIT + 10 glucose strips - English/Italian manual

Unique diagnostic device for self-check of glucose, cholesterol and triglycerides. Multi-parameter technology allows to measure, with the same reflectometer, different parameters by simply changing type of strips: in fact, there is a chip in each pack of strips that is specific for each analyses. Regular measuring of these three parameters is important to **prevent cardiovascular pathologies** if we consider that intolerance to glucose, hyper-cholesterolemia and hyper-triglyceridemia are three of the main cardiovascular risk factors, along with the age, obesity, lifestyle, smoking and arterial hypertension. Multicare is supplied in 3 language carton box (English, Italian, Spanish) with a nylon bag, 10 glucose strips (23950 only), lancet device, 10 lancets, user manual (English, Italian, Spanish as above). Available on request French and Portuguese manual.

23951 MULTICARE DISPLAY KIT
Includes display, multicare kit (23948), glucose strips (23955), cholesterol strips (23957), triglycerides strips (23959), 50 leaflets (English/Spanish).



SELF TEST DEVICE

100% Italian technology

Wide display for easy reading

Multifunction
Each strip box has its own codified chip for safe measurement

GIMA code	STRIPS AND CONTROL SOLUTION FOR MULTICARE
23955	Glucose strips - tube of 25 strips + 1 chip
23956	Glucose strips - tube of 50 strips + 1 chip
23957	Cholesterol strips - tube of 5 strips + 1 chip
23958	Cholesterol strips - tube of 25 strips + 1 chip
23959	Triglycerides strips - tube of 5 strips + 1 chip
23960	Triglycerides strips - tube of 25 strips + 1 chip
23961	Glucose control solution
23962	Cholesterol control solution
23963	Triglycerides control solution



Practical
Device activation by inserting the strip or by pushing the button

Versatile
Blood drop can be applied both on the inserted strip and on the external strip

Innovative
Reduced blood volume 3 microfilters of blood are needed for the test.



TECHNICAL SPECIFICATIONS		
Measuring range:	25-500mg/dl (1.4-28 mmol/l)	Batteries: 2 batteries Cr 2032 3 Volts (32212)
Display:	LCD	Size: 48x104x20.5 mm
Temperature:	10-40°C	Weight: 65 g
Memory:	255 test results with date and time	



STREP-A RAPID TEST, ACCURATE RESULTS IN 5 MINUTES

- 24522 STREP-A - device - box of 20 tests
- 24523 STREP-A - strip - box of 25 tests

Provides accurate results in few minutes, with a simple procedure easy to perform and to interpret. Immediate diagnosis facilitates patient treatment decisions in the same office visit. Room temperature storage. It includes 20 devices or 25 strips tests (single pack), 2 reagents, 20 (25) tubes and swabs. Multilanguage manual and package. For professional use.

Procedure:

- add 4 drops Reagent A, then drops reagent B
- rotate swab 10 times, wait 1 minute, then remove swab while squeezing tube to remove liquid
- add dropper tip to top of tube (24522)
- add 3 drops of specimen (24522)
- immerse strip in specimen (24523)
- read results after 5 minutes



- 24524 STREP-A - twister - box of 25 tests
- Same as code 24522 but with cassette test integrated with the collection tube in one unique device. This make the test procedure easier and quick for the operator. It includes 25 twister cassette, 2 reagents, 25 swabs. Multilanguage manual and package. For professional use.




INFLUENZA A/B TESTS

- 24537 INFLUENZA A/B TESTS - box of 30 tests (3 packs of 10 pcs.)

Influenza A/B is a rapid in vitro test to detect type A and B influenza antigens (nucleoprotein) extracted from the nasal specimens. The test is designed to detect nucleoproteins of Flu-A core subtypes, independent of H and N surface proteins (including H1N1 and H3N2).

A sterile nasal swab allows an easy collection of specimen. Results in 10-15 minutes. Includes rapid influenza A/B test panel, swab, test tube, extraction buffer.

For professional use only.



MONONUCLEOSIS TEST

- 24538 MONONUCLEOSIS TEST - box of 20 tests

A rapid test for the diagnosis of Infectious Mononucleosis (IM) to detect Infectious Mononucleosis heterophile antibodies qualitatively in whole blood, serum or plasma. For professional in vitro diagnostic use only. Results in 5 minutes. It includes 20 cassette tests, 20 droppers, buffers and user manual.



COMBI SCAN 100 URINE ANALYZER AND URINE STRIPS



• **24050 COMBI SCAN 100 URINE ANALYZER** - with printer
 • **33250 PRINTER PAPER 57 mm x Ø 50 mm** - box of 5 rolls
 CombiScan 100 is a small urine test strip reader with modern CCD measurement technology and many features that facilitate daily routine operations.
 To be used with urine strips 24062 and 24063 only.
 Manual in English, German, Italian.
 Software included for PC connection and data transfer.

TECHNICAL SPECIFICATIONS

Measurements: 40 meas./hour in convenient mode 120 meas./hour in fast mode
 Memory: 500 measurements with date and time
 Display: 2 x 24 character LCD-Display
 Printer: Low noise thermal printer with paper roll for more than 300 meas.
 Printout: Current number measurement, Date and time, Serial No. of the meter, patient-I-D, Asterisk on pathological results, Parameters can be printed in SI-units, Conventional units and arbitrary units

Accuracy: Meter adjustment in 5 levels according to individual requirements
 Interfaces: USB and PS2
 Power supply: 100-240 V, 50-60 Hz
 Made in Germany

ANALYTICON COMBI SCREEN® PLUS - urine strips
 Urine strips for visual test or for use with CombiScan 100 (24062/63 only)
 • **24060 COMBI SCREEN® PLUS - 3 parameters** - tube of 100 strips
 - Glucose - Protein - pH
 Results can be read directly on the strip (after about 60 seconds) and compared with the label on the bottle. Supplied in aluminium tube in single carton box. For professional use. Manual in English, French, Greek, Italian, Spanish, Portuguese. Dutch and German available on request. Made in Germany
 • **24061 COMBI SCREEN® PLUS - 10 parameters** - tube of 100 strips
 Same as code 24060 with 10 parameters.
 - Leucocytes - Blood
 - Nitrite - Specific gravity
 - Urobilinogen - Ketones (acetoacetic acid)
 - Protein (albumin) - Bilirubin
 - pH - Glucose
 • **24063 COMBI SCREEN® PLUS - 11 parameters** - tube of 100 strips
 Same as code 24061 plus one parameter: ascorbic acid
 Results can be read on the strips or printed through code 24050.
 • **24062 COMBI SCREEN® SYS - 11 parameters** - tube of 150 strips.
 Same as code 24063 but without ascorbic acid protection PAD on glucose and blood.

OTC, SINGLE PACK URINE STRIPS: UTI, GLU/KET AND 10/11 PARAMETERS

SELF TEST CE 0197

Benefits of OTC Urine Testing

- Convenient: Perform a simple urine test in the comfort of your home, easy to test and read
- Accurate: Obtain reliable results by using the same strips that your doctor uses
- Affordable: Inexpensive testing without the need of laboratory analyzers or home meters
- Fast: Instant result after 2 minutes
- Individually packaged strip in pouch for single testing. This packaging in closed pouch allows 18 months shelf life for strip compared 3 months for remaining strips in an opened canister.
- **24065 GLU/KET STRIP - OTC box containing 3 pouches of 3 strips - CE 0483 - self test**
 Glucose: - Monitor and check excessive glucose in diabetics with screening for a general health check
 Ketone: - Control diabetics by looking for excessive ketones due to low levels of insulin
 - Monitor ketones for sufficient levels during low calorie or low carb diets

- **24066 UTI STRIP - URINARY TRACT INFECTION - OTC box containing 3 pouches of 3 strips - CE 0483 - self test**
 - Effective check Nitrite and Leukocytes in one test for screening of UTI
 - Detect UTI before symptoms are present or before experiencing pain
- **24070 10 PARAMETER STRIP - OTC box containing 3 pouches of 3 strips - CE 0483 - self test**
- **24071 11 PARAMETER STRIP - box containing 5 pouches of 3 strips - CE for professional use**
- **24072 11 PARAMETER STRIP - tube of 50 strips**

These strip are compatible with bayer CLINITEK-100 reader

- Global check/screening for general health condition
- Confirm level of Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose, Ascorbic Acid (24071 and 24072 only) in 1 strip
- Check for healthy nutrition, urinary tract infection, tumors, kidney damage, gout, fever, cirrhosis, viral hepatitis, liver disease, congestive heart failure, jaundice, diabetes and acute renal failure.



WARNING: SHORT EXPIRY DATE FOR TESTS. Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

CARDIAC MARKER

• 24525 CARDIAC MARKER - 3 parameters:

CK-MB, Troponin, Myoglobin - box of 20 tests
It allows to check cardiac status and stroke events in just 15 minutes.

Ideal to keep under control patients affected by cardiopathies or after heart attack or surgery. Acute Myocardial Infarction-MI (Heart Attack) is responsible for 1 in 5 deaths in the U.S. After 6 hours of AMI onset, the damage to the cardiac muscle is irreversible and permanent.


This illustrates that early and accurate diagnosis of AMI is essential. ECG is the most widely used method of MI diagnosis. Unfortunately, as many as 25% of AMI patients have a normal ECG. GIMA offers an array tests for Biochemical and Protein markers of MI, such as Creatine Kinase Isoenzyme MB (CK-MB), Myoglobin (MB) and Cardiac Troponin I (cTnI) in combination. These parameters are well recognized and accepted worldwide. It includes 20 droppers, 20 cassettes and 1 bottle with dilution solution. For professional use.



D-DIMER TEST

• 24526 D-DIMER TEST - box of 20 tests

The presence of D-Dimer in blood samples is an indicator of various coagulation disorders, including deep venous thrombosis (DVT), pulmonary embolism (PE) and atherosclerosis. D-Dimer assay is a widely used and simple exclusion method of DVT and PE. This test allows an easy prevention: a negative result indicates a good health status, while high concentration indicates increased activity of fibrin and blood coagulation. The sensitivity of the Ameritek D-Dimer test is 97% (174/178) and the specificity is 92% (114/122) when compared to the golden standard Dade Behring Stratus CS - DDMR Immunoassay System. It includes 20 droppers, 20 cassettes and 1 bottle with dilution solution. For professional use.



MALARIA TEST


• 24547 MALARIA TEST - box of 40 tests

The Malaria P.f. Rapid Test Device (Whole Blood) is a rapid chromatographic immunoassay test for the qualitative detection of circulating plasmodium falciparum in whole blood. It has been compared with traditional thick or thin blood smears microscopic analysis, the correlation between the two systems is over 99,0%. Package includes test devices, buffer, disposable specimen dropper and multilanguage manual, lancets are not included.

For professional use.



24547



SEMI-QUANTITATIVE CRP TEST

• 24527 CRP TEST - box of 20 tests

Measurement of CRP in whole blood, plasma or serum has two purposes:

- 1) hs-CRP levels help in assess the risk of cardiovascular and peripheral vascular disease. Higher hs-CRP values have highest risk of cardiovascular disease.
 - 2) CRP levels help identify and keep track of infections and diseases that cause inflammation (lymphonoma, lupus, giant cell arteritis, osteomyelitis, rheumatoid arthritis,).
- Use fingerstick whole blood or serum/plasma.
 - Results with in 5 to 10 minutes.
 - Strong intensity colour bands.
 - Built in performance control and reference band
 - Room temperature storage.

It includes 20 cassette tests with droppers. For professional use.




24527

FOB-FECAL OCCULT BLOOD TEST FOR EARLY DETECTION

• 24530 FOB TEST - box of 25 tests

Provides accurate results in few minutes. It facilitates patient treatment decision early. Simple procedure, easy to perform and interpret. No additional equipment required. It includes 25 pcs of collection stick, device test and one buffer bottle.

For professional use



24530

ALLERGY SCREENING TESTS: CAT, DOG AND DUST



24580

- 24580 CAT ALLERGY TEST - box of 1 test
- 24581 DOG ALLERGY TEST - box of 1 test
- 24582 DUST ALLERGY TEST - box of 1 test

SELF TEST

Simple, reliable and clinically proven screening test for allergies:

- Result in 15 minutes
- Excellent correlation to lab test
- Easy to use a fingerstick whole blood
- For self testing
- Each test includes: test device, vial of test solution, fingerstick device, user manual




MICROALBUMIN URINE SCREENING TEST - KIDNEY DISEASE

• 24539 MICROALBUMIN URINE TEST - box of 1 test

An easy-to-use visual test to detect elevated levels of albumin in urine, an early warning sign of kidney disease. Early detection is key: two main causes of kidney disease are **diabetes & high blood pressure**.

- Excellent correlation to lab test
- Results in 3 minutes
- Simple, easy to use urine test
- Sensitive; detects 18 mg/L of albumin in urine (same as professional tests)
- Contains everything needed for consumer to perform test
- CE approved as self-testing device
- American Diabetes Association and American Kidney Foundation recommend that those with diabetes test annually.



SELF TEST

COLON DISEASE TEST, FOR SELF TESTING AT HOME

SELF TEST

• 24531 COLON DISEASE SELF TEST

Simple home test for detecting the early warning signs of colorectal disease, even before symptoms appear.


Contains 5 test pads, 1 positive control package, results card and 1 instruction manual.

- Easy throw in the bowl test for hidden blood in the stool
- No laboratory processing required, results in 2 minutes
- The American Cancer Society recommends annual testing after age of 50
- With early detection and treatment colorectal cancer is curable

CE 0197
FDA Approved



24531



Step1: after a bowel movement drop a test pad into the toilet bowl.

Step2: observe the test pad. If blood is present a blue colour will appear within two minutes

HELICOBACTER PYLORI TEST FOR ULCER SCREENING

CE

• 24528 HELICOBACTER PYLORI TEST - box of 40 tests

Helicobacter Pylori is the leading cause of gastric ulcers, this easy to use test allows a quick diagnosis in about 10 minutes. Easy to perform and interpret, no additional equipment required. For professional use.

Diagnosis of H. pylori infection in adults over 18 years of age.

• 24529 HELICOBACTER PYLORI SELF TEST - box of 1 test

Same as 24528 but packed in single box including one automatic lancet, single buffer and user manual.

This product is provided with self testing certification to be sold to end users in OTC market.

SELF TEST



24529



24528



PROSTATE SCREENING - PSA ONE-STEP TEST

• 24533 PSA TEST - box of 40 tests - CE - for doctor/hospital use

Provides accurate results in few minutes.

- Over 98% correlation to lab tests
- Simple, easy to use whole blood test
- Detects PSA as low as 4ng/ml
- Contains everything needed to perform test
- For professional use

• 24532 PSA SELF TEST - box of 1 test - CE 0197 - for self test

Same as 24533 but packed in single box including one automatic lancet, single buffer and user manual. This product is provided of self testing certification to be sold to end user.

SELF TEST



24532



24532

VISUAL STRIPS FOR SELF TEST GLUCOSE AND CHOLESTEROL

• 24520 GLUCOSE STRIPS - box of 30 pcs.

• 24521 CHOLESTEROL STRIPS - box of 2 pcs.

A new, inexpensive, CE marked line of strips for visual reading of cholesterol and glucose in whole blood. **No meter is needed to measure the results.** The special color chart allows an easier comparison of the colors for a better quantification of the result. Cholesterol is offered in a 2 tests tube package inclusive of fingertips, color chart and instructions for use.

Glucose is offered in a 30 tests/tube package inclusive of colour chart and instructions for use. Box and instructions in English, German, Italian.

CE 0197



24521



24520

SELF TEST

BREATH ALCOHOL TESTER

24500 ALCOHOL TESTER

Compact device to measure your alcohol concentration in your breath. Digital/graphic interface and animation, one-hand operation. % BAC and mg/l values.

Supplied in carton box with 3 blowing tubes Ø 7 mm and one cleaning cap.


- 24501 BLOWING TUBE - box of 12 pcs.
- 32200 "AAA" BATTERIES - blister of 4 pieces
- 24503 COUNTER DISPLAY - Italian
- 24504 COUNTER DISPLAY - English



Large display to read exact value of your breath (3 colour backlight)

SELF TEST

Counter display fits 6 tester



24500

TECHNICAL SPECIFICATIONS

Range: Low & 0.15 ~ 1.50 % BAC or 0.075 ~ 0.75 mg/l
 Accuracy: ± 0.2% BAC full range under alert level
 Battery: 3 x 1.5V AAA batteries (not included)
 Display: Graphic & Digital display with 3 decimal points (eg. 0.025%BAC)
 Tolerance compensation: Pre-alert function together with tri-colour backlighting
 Breath sampling: less than 5 seconds continuous breath sample
 Recycle time: about 20 seconds
 Unit Dimensions: 107 x 45 X 25 mm approx.
 Weight: 50g (without batteries)
 Alert level calibrated: 0.5 % BAC or 0.25mg/l



PROFESSIONAL ALCOHOL TESTER

- 24491 FIT-3 ALCOHOL TESTER with printer
- 24494 BLOWING TUBES - box of 100 pcs.
- 24495 PAPER ROLLS for 24491 - box of 10
- 24496 PRINTER CARTRIDGE for 24491

Professional fuel cell sensor breath alcohol tester with printer and touch screen. The fuel cell is a porous disk with a thin layer of platinum black on both faces and any alcohol that comes in contact with the cell's surface is broken down chemically. For each molecule of alcohol that is broken down a given number of electrons are freed during the process. The numbers of electrons generated are in proportion to the amount of alcohol.

Alcohol tester¹ is supplied with integrated printer, 25 mouth pieces², 1 large mouth piece³, power supply⁴, 12V power supply⁵, wireless module⁶, USB car adaptor⁷, wireless transmission module⁹, carrying pouch¹⁰, 4 rolls of paper¹¹, charger¹².



Size: 188x78x33 mm
Weight: 300 g

TECHNICAL SPECIFICATIONS

Display:	Colour 3.2 TFT LCD with touch screen and graphic waveform
Keyboard:	5 keys operation: "Up", "Down", "Left", "Right", "OK"
Power Supply:	Rechargeable AA battery (about 2500 mAh)
Heat warm up:	Chamber & Sensor warm up automatically
Auto flow monitoring:	Flow detect automatically
Storage environment:	-40 Degree C to +85 Degree C.
Operating environment:	-20° to +70°. (optimal -20° to +40°)
Record PC download:	PC Download through RS-232 cable min. 10,000-18,000 memories with full information: drive name, vehicle number, inspector name, test result with stamp data and time, location, calibration data, ...
Memories:	
Measurement precision:	Standard deviation +/-0,025 mg/L or +/-5% of the measurement Drift: typically <0,8% of the measurement value/month

POCKET ALCOHOL TESTER

- 24499 POCKET ALCOHOL TESTER

Pocket, easy to use alcohol tester

TECHNICAL SPECIFICATIONS

Range: Low & 0.2 ~ 1.5 % BAC
 Accuracy: ± 0.2% BAC at the calibrated point (alert level)
 Power input: 2 x 1.5V AAA batteries (not included)
 Display: Digital display
 Tolerance compensation: pre-alert function when close to alcohol limit

Auto power off: 2 minutes
 Breath sampling with auto flow detect: about 6 seconds continuous breath sample
 Cycle time: about 20 seconds
 Alert level calibrated: 0.5% BAC



POCKET, ONLY
8.5x3x1 cm




SALIVA ALCOHOL STRIPS

- 24502 SALIVA ALCOHOL TEST STRIPS - box of 25 strips

GIMA Saliva Alcohol Test Strip is a rapid highly sensitive test for the determination of relative blood alcohol concentration in saliva specimens.

- Results in just 2 minutes
- Color chart provides approximate relative blood alcohol concentrations (BAC) from 0.02% to 0.30%
- No instrumentation is required
- Saliva testing avoids urine collection
- CE Marked for sale in the European Union
- Short expiry date (average 8-10 months)



SELF TEST

DRUG TEST FOR HOME CARE USE: URINE OR "WIPE" TEST

- 24560 MULTIDRUG SELF TEST urine

Multi drug test for home testing, supplied with urine container and manual.

Easy to use, provides accurate results in few minutes for the following drugs: amphetamine(AMP), cocaine(COC), marijuana cannabis (THC), methamphetamine (METH) including ecstasy, opiates (MOR), including heroine, phencyclidine (PCP).

Supplied Italian/English package and manual. Test can be offered in any language (box+manual) minimum 1000 pcs.

SELF TEST



24560

- 24564 DRUG KOP - box of 1 test
- 24565 DRUG KOP - box of 10 tests

A new colorimetric test to screen out quickly and accurately the major substances contained in illicit drugs and some pharmaceuticals. This test measures and detects the following 21 DOA with concentration of 5 mg: Amphetamine, Benzphetamine, Cocaine, Codeine, Dextromethorphan, Diphenhydramine, Doxylamine, Heroin, Hydrocodone, Ketamine, LSD (containing a secondary amine), Mephentermine, Methadone, Methamphetamine, Morphine, Opium, Phencyclidine (PCP), Promethazine, Thebaine, Tryptamine, Tyramine.

Test can be performed in surfaces and drinks, it does not require reagents and gives immediate results. Very stable from - 15°C to 50°C.



SELF TEST

DRUG TEST READER SYSTEM - FOR URINE DRUG CUPS/PANELS

- 24570 DRUG READER - Italian software
 - 24571 DRUG READER - English software
 - 24574 MULTI DRUG PANEL - 8 drugs
 - 24576 MULTI DRUG CUP - 7 drugs + S.V.T.
 - 24577 MULTI DRUG CUP - 8 drugs + S.V.T.
- Simply scan one of the three multi drug test on drug reader and get a complete analysis at your PC.
- automated analysis of "visual read" test panels
 - automatic test menu recognition (2D barcode)
 - in-process system quality control check
 - image capture, data management/export functionality, archived results
 - validated for use with multiple formats/test menus
- Main benefits:
- automation - eliminates results interpretation variability/subjectivity
 - accurate - performance advantages vs. visual interpretation
 - speed - scanning process requires 15-30 sec.
 - flexible - multiple test formats/menu, semi-portable, simple installation, Windows compatible.

DRUG	24574	24576	24577	Detection time
MTD	•	•	•	5-7 days
MDMA	•	•	•	5-7 days
MET	•	•	•	3-5 days
COC	•	•	•	2-5 days
THC	•	•	•	-
AMP	•	•	•	2-5 days
MOP	•	•	•	5-7 days
BUP	•	•	•	5-7 days
S.V.T.				
OX		•	•	
S.G.		•	•	
PH		•	•	
NiT		•	•	
Glut		•	•	
CRE		•	•	



MULTI DRUG "CUP" TEST

- 24546 MULTI DRUG "CUP" TEST - urine box of 25 tests for 7 drugs

Detects following parameters: (MTD/MDMA)+<COC>+<AMP500>+<MET500>+<THC>+<MOP>. For professional use.

- Easy to use - no urine handling, just insert the key into the active test and read results in 5 min.
- Tightly sealed urine cup - no leaking or spilling of urine specimen
- Unique design - ensures that urine used in the test, will be separated from residual urine, that will remain uncontaminated, to be used for a further test when necessary
- Key packed separately, ensures security that only the technician can begin the test

Test panel is tightly fastened to the cup, protecting it from contamination



SALIVA DRUG TEST

- 24555 MULTI DRUG SALIVA TEST oral fluid (saliva) - box of 25 tests for 6 drugs

Multi drug tests provide accurate results in 10 min. for the following drugs: amphetamine, cocaine, marijuana cannabis, methamphetamine, opiates, phencyclidine. Cut-off levels meet SAMSHA standards. For hospital and professional use. Multilanguage manual: English, Italian, French, German, Spanish, Portuguese.



Collect oral fluid with sponge

CE
FDA Approved

SCREENING TEST CUP - 24546		
DRUG	CUT OFF NG/ML	DETECTION TIME*
MOP	300	5-7 days
COC	300	2-5 days
THC	50	-
AMP	500	2-5 days
MDMA	500	5-7 days
MTD	300	5-7 days
MET	500	3-5 days

ORAL FLUID SCREENING TEST DEVICE - 24555		
DRUG	CUT OFF NG/ML	DETECTION TIME
COC	20	10 min - 24 h
THC	12	up to 14 h
AMP	50	10 min - 72 h
MET	50	10 min - 72 h
PCP	10	-
OPI	40	1 h - more days


* Drug use over detection time not measured

DRUGS OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

CE FDA Approved

- 24540 COCAINE TEST - strip - urine - cut-off: 300 ng/ml
 - 24541 MARIJUANA TEST - strip - urine - cut-off: 50 ng/ml
 - 24542 OPIATES TEST - strip - urine - cut-off: 300 ng/ml
- Code 24540-41-42 are sold in box of 50 strips individually packed in single pouch.

- 24550 MULTI DRUG PANEL TESTS urine box of 25 tests for 6 drugs
 - 24551 MULTI DRUG PANEL TESTS urine box of 25 tests for 12 drugs
- Multi drug tests provide accurate results in 10 min. on the following drugs: amphetamine*, cocaine*, marijuana cannabis*, methamphetamine*, opiates*, phencyclidine*, barbiturates, benzodiazepines, ecstasy, morphine300, methadone, tricyclic antidepressants. *Code 24550 contains tests for drugs with asterisk. Cut-off levels meets SAMSHA standards. All products are supplied with multilingual manual: English, Italian, French, Spanish, German, Portuguese.



SCREENING TEST CUP			
Drug	Cut off ng/ml	Detection time	
MOP	24550 300	24551 300	5-7 days
COC	300	50	2-5 days
THC	50	50	5-8 days
AMP	1000	1000	2-5 days
MDMA	-	500	5-7 days
MTD	-	300	5-7 days
MET	1000	1000	3-5 days
BZO	1000	300	7-10 days
BAR	-	300	2 days - 4 weeks
PCP	25	25	-
OPI	2000	2000	5-7 days
TCA	-	1000	-

PROTECTIVE GOGGLES



- 25661 POLYSAFE MEDICAL - bulk
 - 25660 POLYSAFE MEDICAL - with box
- They are a vital part of any infectious disease control program for eye protection. They provide a wide field of vision and excellent optical quality. Designed with top and side shields to protect eyes from every angle. They can fit all head sizes comfortably, and have superior optical clarity. Autoclavable. Norms EN 166.



- 25659 GIMASAFE - box of 10 pcs
- Transparent protective goggle.

VISOR HEAD MASK


- 25656 VISOR HEAD MASK
- Great features**
- faceshield with 90° max rotation
 - inner and outer surface treated against mist
 - superior coverage to protect against droplets and splashes of liquids
- Superior materials**
- elastic retention band in anallergic material
 - rotation systems in anodized aluminium
 - soft support in anallergic foam
 - faceshield in heat-formed acetate with anti-fog treatment
- Easy and comfortable use**
- fitting headbands adjustable by velcro fasteners
 - usable over ophthalmic glasses



- 25644 SAN DIEGO GOGGLES - black - fog resistant
 - 25645 SAN DIEGO GOGGLES - blue - anti scratch
- Temple lenght and inclination adjustment by snap motion. Protection against high speed particles - low energy and lateral protection. Protection against optical radiation to standard EN 170.



- 25664 MILLENIA GOGGLES - black - fog resistant
- The best quality of protection goggle according to a positionable and innovative design.



- 25663 HIGH PROTECTION GOGGLES - fog resistant and antiscratch
- Protection against: droplets and splashes of liquids and high-speed particles at extreme temperature. Lateral protection. Protection against large dust particles (>0.5 μ). Protection against optical radiation to standard EN 170. Resistance to surface damage by fine particles. Resistance to fogging of oculars.



- 25646 GOGGLES KIT - fog resistant
 - 25647 SPARE SHIELDS - box of 20
- Pratical kit composed of 5 supports in different colours and 20 antifog shields. Usable over ophthalmic glasses.

VINYL AND LATEX GLOVES

VINYL GLOVES: Sensitivity and protection.

Ideal gloves for hospital attendants, for everyday use, they warrant protection and sensitivity. Medium weight, designed for many applications, including laboratory procedures and hospital ward use. Packed in dispenser box 100 units. EN455-EN347-AQL<1.5



With Powder	Powder free	Powder free elasticized	VINYL GLOVES - NON STERILE (box of 100)
25701	25692	25716	Small size - NON sterile
25702	25693	25717	Medium size - NON sterile
25703	25694	25718	Large size - NON sterile

LATEX GLOVES ambidextrous - with hem
Pure latex gloves, non sterile, hypoallergenic, packed in dispenser box 100 units.

EN455-EN347-AQL<1.5



With Powder	Powder free	LATEX GLOVES NON STERILE (box of 100)
25706	25711	Small size - NON sterile
25707	25712	Medium size - NON sterile
25708	25713	Large size - NON sterile
With Powder	Powder free	LATEX GLOVES EXTRA LONG NON STERILE (box of 50)
-	25691	Medium extra long size 400 mm

STERILE SURGICAL GLOVES

Latex sterile surgical gloves with low content of latex proteins. The micro textured surface assures a very good grip in any different condition of use. Lightly powdered. Anatomic shape assuring a very good comfort. Sterilized by ETO.



With powder	Powder free	LATEX SURGICAL GLOVES - STERILE box of 50 pairs
25685	25695	Size 6.5 - sterile
25686	25696	Size 7 - sterile
25687	25697	Size 7.5 - sterile
25688	25698	Size 8 - sterile
25689	25699	Size 8.5 - sterile



DRACULA



For minimum 500 pcs. we can put your logo here

• **25700 DRACULA**
Universal holder for gloves and tissue boxes, with little sharp teeth automatically adjust to tightly hold any size of box. It's an universal holder that allows one-hand operation.

LATEX FINGER COTS

• **25941 LATEX FINGER COTS - medium**
• **25942 LATEX FINGER COTS - large**
Latex finger cots to unroll 1 finger. No sterile, single use, for fingers protection, with talc, no antistatic.
minimum order: 1000 pcs. (10 bags of 100 pcs.)



NYTRILE GLOVES - NO POWDER

Disposable nitrile examination gloves, powder free and latex free, conform EN455 - EN374 - AQL<1,5 Suitable for end users with sensible skin, medical, labs, final industry use. They assure a great protection. Veloforte model is particularly sensible and elastic. Light blue.



Great protection

Standard	Veloforte	NYTRILE GLOVES - NON STERILE (box of 100)
25680	25679	Small size - non sterile
25681	25683	Medium size - non sterile
25682	25684	Large size - non sterile



COPOLYMER GLOVES

COPOLYMER GLOVES ON PAPER
Examination gloves, elasticized, medium size, for gynecology use.
• **25714 COPOLYMER GLOVES - non sterile**
minimum order: 10 pack of 100 pieces = 1.000 pcs.
• **25709 COPOLYMER GLOVES - sterile**
Sterilized by ETO. AQL <1.5
minimum order: pack of 100 pieces



COTTON GLOVES

White	Green	Size	Length cm		Weight g	
			min	max	min	max
25864	25874	Size 6.5	20	22	18	20
25865	25875	Size 7	21	23	19	21
25866	25876	Size 7.5	23	25	21	23
25867	25877	Size 8	24	26	23	25
25868	25878	Size 8.5	25	27	25	27

Minimum order: 1 box of 10 pairs



VISOR SHIELD "PROTECTOR"

• **25642 VISOR SHIELD "PROTECTOR"**
The ultra-light face shield that combines function with style:
protection: a waterproof plastic barrier will protect you
comfort: fast and easy to wear and remove
convenience: reusable and replaceable shield
optical quality: superior optical quality, the best you can buy
Medium adjustable size: fits over 90% of healthcare professionals.
Fits over glasses. White colour.
Supplied with 2 face transparent shields + 1 head shield.
• **25643 REPLACEMENT SHIELD KIT**
Includes 2 face transparent shields + 1 head shield



Quality product: made in USA



25642



99% FILTERING MASKS 3 PLY - WITH LACETS OR LOOPS

99% FILTERING MASKS

- Lightweight and comfortable for soft skin
- Size 17x9 cm (adult) or 12x9 cm (pediatric)
- Conformable nose band for anatomical fitting
- Fiber glass free, filtration efficiency of BFE +99% at 1.0 micron
- Antifog, latex free
- Suitable for hospital, dental and ambulatory use
- With lacets or ear loops
- EN 14683 conform

8 languages box



25658

Mask with ear loop



25657

Mask with lacets



Masks with ear loops

GIMA code	99% FILTERING, 3 PLIES MASKS box of 1,000 (20 dispenser of 50 pcs.)
25653	Adult 3 plies masks - ear loop - green
25658	Adult 3 plies masks - ear loop - light blue
25673	Adult 3 plies masks - ear loop - pink
25657	Adult 3 plies masks - lacets - green
25674	Pediatric 3 plies masks - ear loop - white

99% FILTERING MASKS MADE IN ITALY - EN 14683:2005



CE
93/42

• 25648 MASKELLA - face mask with ear loops - light blue

- 3 layers of non-allergenic woven-non woven
- Size 90x180 mm
- Weight 3 g
- Tensile strenght extra-long ties: 80 cm
- Comfortable, shapeable nose clip
- Filtering efficiency > 99%
- Resistance to breath (MM H₂O/cm²) < 2
- Does not contain glass fibres
- Does not contain latex
- Non sterile

min. order : 18 packs of 50=900



CE
93/42

• 25649 MASKERITA - surgical face mask

- 3 layers of non-allergenic woven-non woven
- Size 90x180 mm
- Weight 3 g
- Tensile strenght extra-long ties: 80 cm
- Comfortable, shapeable nose clip
- Filtering efficiency > 99%
- Resistance to breath (MM H₂O/cm²) < 2
- Does not contain glass fibres
- Does not contain latex
- Non sterile

min. order: 12 packs of 50=600



CE
93/42

• 25652 MASKOP

- Surgical mask with wrap-around visor
- Anti-misting, antistatic wrap-around visor
- 4 layers of non-allergenic woven-non w.
- Weight 10.5 g
- Atraumatic and conformable nose-clip
- Filtering efficiency > 99%
- Resistance to respiratory flow (MM H₂O/cm²) < 2.5
- Latex and glass fibres free. Not sterile
- Fluid permeability: impermeable
- External/Internal polyethylene anti-misting strip

min. order: box of 50

2 PLY MASKS

25651



• 25651 2 PLY MASKS - light blue box of 2,000 pcs. (40 dispenser of 50 pcs)

- Made in soft TNT
- Size 17x9 cm with ear loops
- Comfortable nose band
- Antifog, latex free
- EN 14683 conform

3 PLY MASKS 95%

25650



• 25650 3 PLY MASKS - green box of 1,000 pcs. (20 dispenser of 50 pcs)

- 95% filtering mask
- Lightweight and comfortable for soft skin
- Size 17x9 cm with ear loops
- Conformable nose band for anatomical fitting
- Fiber glass free, filtration efficiency of BFE +95% at 1.0 micron
- Antifog, latex free
- EN 14683 conform

MOULDED MASKS



• 25677 MOULDED CONE MASKS - light blue dispenser of 50 pcs.

- 100% latex free, including the elastic headband
- Splash repellent barrier
- Flared edge prevents irritation
- Soft, flexible, aluminum nosepiece adjusts for custom fit
- Soft, odorless, fiberglass free
- High breathability
- BFE exceeds 99% at 4.5 microns

FOLD-FLAT RESPIRATORS FFP1, FFP2, FFP3



Facial filtering respirators to protect against solid (non-oil) and liquid (oil) volatile particles. With exhalation Valve (excluding 25665). Double elastic with buttonhole allows to hang the mask to the neck when it's not worn. Aluminium adjustable nose piece soldered with high resistance and internal Polyethylene (PE) Foam nose pad. Single packed guarantees total hygiene and

easy maintenance. The shape of the respirator has been developed to reduce difficulties of respiration and to improve efficiency in the dialogue. Innovative design allows a perfect view. Wide adaptability to various face contours. Wide filtering area to allow proper ventilation. Light weight. Conforms to CE standard - EN 149: 2001

GIMA code	PARTICULATE FOLDABLE RESPIRATORS	minimum order
25666	Particulate respirator - FFP1 - with valve	box of 20
25665	Particulate respirator - FFP2	box of 20
25670	Particulate respirator - FFP2 - with valve	box of 10
25671	Particulate respirator - FFP3 - with valve	box of 10

	TECHNICAL SPECIFICATIONS			
	25666	25665	25670	25671
Protection:	FFP1	FFP2	FFP2	FFP3
Minimum filter efficiency	78%	92%	92%	98%
Exposure limit fine particles up to	4xTLV	12xTLV	12xTLV	50xTLV
Valve	yes	no	yes	yes
Vertically foldable	yes	yes	yes	yes

MATERIALS	
Filter: Polypropylene Fabric	Elastics: latex-free
Valve: Plastic rigide (excluding 25665)	
Valve Membrane: Synthetic Rubber (excluding 25665)	
Nose Piece: Aluminium adjustable	Nose Pad: Polyethylene (PE) Foam

NORMS: CEI EN 149: 2001
AC 538100 05-06-2008

BOUFFANT, SURGEON CAP

- 25675 BOUFFANT CAP - light blue
Minimum order: 1,000 pcs (10 bags of 100 pcs)
- made of light weight PPSB
- breathable, soft, comfortable
- elastic headband
- medium size 21"



25675


- 25678 SURGEON CAP green
Minimum order: 1,000 pcs (10 bags of 100 pcs)



25678

SHOE COVERS


- 25676 SHOE COVERS - light blue
Minimum order: 2,000 pcs (20 bags of 100 pcs)
- made of PPSB, breathable
- comfortable
- with elastic



25676

BIBS

- 25672 BIBS - roll of 80 - 50x60 cm
Made of 2 plies, one of highly absorbent cellulose tissue and one of waterproof polyethylene. Extremely soft and comfortable, they assure an excellent protection to the patient. Each roll has 80 bibs with pre-cut and comfortable ribbons or bi-adhesive.



Minimum order: box of 6 rolls

FOLDED NAPKINS

FOLDED NAPKINS
Pre-coated folded napkins, made of 2 plies of pure wadding of cellulose partially laminated on 1 ply of polyethylene. The special features of materials offer high absorbency and waterproofness, to assure a high comfort to the patient.




GIMA code	FOLDED NAPKINS - PACK OF 500 PCS.	
25230	Napkins 33x45 cm	- green colour ■
25231	Napkins 33x45 cm	- light blue colour ■
25232	Napkins 33x45 cm	- pink colour ■
25233	Napkins 33x45 cm	- white colour
25240	Plastic dispenser	
25239	Chain for napkins	

OTHER COLOURS AND SIZES AVAILABLE FOR QUANTITIES

SURGICAL CAPS



- 26166 GREEN CAP
Other colour on request. Minimum order 1000 pieces.
- 26165 FANTASY CAP - light blue
- 26167 FANTASY CAP - beige
- 26168 FANTASY CAP - red
100% cotton surgical cap with fantasy pictures or green.



DOCTOR'S WHITE AND SURGICAL COATS, JACKETS, TROUSERS



COATS, TROUSERS AND JACKETS
Professional and stylish line with all sizes for men and women, made in high quality 100% cotton.
Available with buttons or press studs.
Conform to 89/686/CEE.

DOCTOR'S WHITE COATS 100% COTTON						JACKETS				
Size			Woman	Man	Unisex with press stud	Size	White	Green	Light Blue	Fantasy
Europe	USA	France								
40	XS	36	26110	-	26080	XS	26124	26126	26128	
42	S	38	26111	-	26081	S	26130	26140	26150	26170
44	S	40	26112	26098	26082	M	26131	26141	26151	26171
46	M	42	26113	26101	26083	L	26132	26142	26152	26172
48	M	44	26114	26102	26084	XL	26133	26143	26153	26173
50	L	46	26115	26103	26085	XXL	26134	26144	26154	-
52	L	48	26116	26104	26086	TROUSERS				
54	XL	50	26117	26105	26087	XS	26125	26127	26129	-
56	XL	52	-	26106	26088	S	26135	26145	26155	-
58	XXL	54	-	26107	26089	M	26136	26146	26156	-
						L	26137	26147	26157	-
						XL	26138	26148	26158	-
						XXL	26139	26149	26159	-

THEATRE GOWNS - GREEN	
Size	Man-woman
52-56	26162
58-62	26163

JACKETS WITH PRESS STUD AND TROUSERS

JACKETS AND TROUSERS



JACKET WITH STUD AND TROUSERS
Made in 100% pure drill cotton jackets are offered in white colour from stock or in 6 different colours in 30 days.

JACKETS WITH PRESSURE STUD						
Size	Woman			Unisex		
	White	Green	Violet	White	Blue	Light blue
XS	21101	21111	21121	21151	21161	21171
S	21102	21112	21122	21152	21162	21172
M	21103	21113	21123	21153	21163	21173
L	21104	21114	21124	21154	21164	21174
XL	21105	21115	21125	21155	21165	21175
XXL	21106	21116	21126	21156	21166	21176
Size	TROUSERS - UNISEX					
	White	Green	Violet	White	Blue	light blue
XS	26125	21211	21221	26125	21261	21271
S	26135	21212	21222	26135	21262	21272
M	26136	21213	21223	26136	21263	21273
L	26137	21214	21224	26137	21264	21274
XL	26138	21215	21225	26138	21265	21275
XXL	26139	21216	21226	26139	21266	21276

DENTAL JACKETS

DENTAL JACKETS



DENTAL JACKETS
Short sleeves round-necked jacket with 2 buttons on right shoulder and 5 hidden buttons. One pocket and one breast pocket. Green colour 100% pure cotton

DENTAL JACKET	
Size	Unisex green
XS	21000
S	21001
M	21002
L	21003
XL	21004
XXL	21005

ANTISTATIC CLOGS - AUTOCLAVABLE AT 135°C

GIMA PROFESSIONAL CLOGS

- active arch support to stimulate circulation
 - antistatic material directly injected in the polymeric blend
 - perforated sole for breathability and transpiration
 - antislip grooved sole safe on wet floor
 - machine washable, sterilizable at 134°C
 - available also with strap for extra security (26471-26482)
 - individually packed in carton box
 - all models (with or without pores, with or without straps) are available in any size or colour (including red).
- Minimum order for products not coded: 100 pairs.



26331 - 26342

UPPER WITHOUT PORES

UPPER WITH PORES

26451 - 26462

26231 - 26242

WITH MOVABLE STRAPS WITH OR WITHOUT PORES



EN ISO 20347 A-E

Made in Italy

AVAILABLE IN OTHER COLOURS/SIZES (minimum order 100 pairs)

Size	UPPER WITHOUT PORES						UPPER WITH PORES				
	White	Green	Electric Blue	Orange	Fuchsia	White with strap	White	Yellow	Blue	Violet with strap	Green with strap
34-35	26311	26331	-	26371	-	26231	26411	-	26451	-	26471
35-36	26312	26332	-	26372	-	26232	26412	-	26452	-	26472
36-37	26313	26333	-	26373	-	26233	26413	-	26453	-	26473
37-38	26314	26334	26354	26374	26394	26234	26414	26434	26454	26254	26474
38-39	26315	26335	26355	26375	-	26235	26415	-	26455	26255	26475
39-40	26316	26336	26356	26376	-	26236	26416	-	26456	26256	26476
40-41	26317	26337	26357	26377	-	26237	26417	-	26457	26257	26477
41-42	26318	26338	26358	26378	26398	26238	26418	26438	26458	26258	26478
42-43	26319	26339	26359	26379	26399	26239	26419	26439	26459	26259	26479
43-44	26320	26340	26360	26380	26400	26240	26420	26440	26460	26260	26480
45-46	26321	26341	-	26381	-	26241	26421	-	26461	-	26481
47-48	26322	26342	-	26382	-	26242	26422	-	26462	-	26482



DISPOSABLE GOWNS POLYTHENE GOWNS SURGICAL GOWNS

DISPOSABLE SURGICAL GOWNS - green

Disposable non woven gowns, suitable when comfort and light protection are required. Colour: green Grams: 25 g/sqm One size: 110x136 cm With cotton cuffs (26100 only)



26100

DISPOSABLE POLYTHENE GOWNS - green

Surgical polythene gowns with neck and waist laces and cotton cuff. Made in waterproof TNT. One size 115x140 cm Available sterile or not sterile. Made in Italy.



DISPOSABLE SURGICAL SPECIALISTIC GOWNS - light blue - sterile

Specialistic, water-proof, pre-folded surgical gowns for easy, programmed dressing: - cotton cuff - double, internal seam - neck velcro closure - 4 waist straps for internal and external closure Includes hand wipe Made in Italy.



GIMA code	DISPOSABLE GOWNS	Minimum order	GIMA code	POLYTHENE GOWNS	Minimum order	GIMA code	SURGICAL GOWNS STERILE	Minimum order
26099	Disposable Gowns - sterile*	box 50	26092	Poly Gowns - sterile	box 50	26095	Size L - 130x150 cm	box 20
26100	Disposable Gowns - non sterile	bag 100	26093	Poly Gowns - non sterile	box 100	26096	Size XL - 140x160 cm	box 20

*Sterilized by ETO

GIMA SECURITY SYSTEM TO FUSE NEEDLES UP TO Ø 1.6 MM

• 25820 SECURITY SYSTEM - needle destroyer

To avoid accidents from used injection needles, GIMA has developed and patented Security System. This Needle Destroyer generates a temperature of 1,000°C, hot enough to destroy stainless steel and all known bacteria/ viruses. It fuses in 1-2 seconds needles up to a diameter of 1.2 mm and with some more seconds needles up to 1.6 mm. Safe for nurses, doctors, cleaning staff and patients. Most of accidents by sticks with needles take place in the following circumstances that are avoided using the Needle Destroyer: placing the needle in its protective sheath, discarding a needle from a used syringe, carrying hospital residues. Destroyed needle material is sterilized and stored in a drawer at the bottom of the unit. Long life and easy maintenance. It is sturdy, over 4 kg weight, and manufactured in best quality material. Only necessary maintenance is to change fusing spare parts every 4,000 needles.

TECHNICAL SPECIFICATIONS		
Size: cm 22x14x h 9	Weight: 4 kg	Power: 400 Watt Fuses: 2x2,5 A T
Operating voltage: 220 V - 50 Hz	Fusion temperature: 1,000°C	Norms: CEI 62-5 - IEC 601-1
Made in Italy: GIMA S.p.A.		





• 25821 FUSING SPARE PARTS minimum order: 2 pcs.

SECURITY SYSTEM ELIMINATES ACCIDENTAL CONTAGION OF INFECTIOUS DISEASES CAUSED BY NEEDLE-STICKS

TOURNIQUETS IN LATEX, RUBBER-COTTON AND LATEX-FREE



- 25722 JETPULL 2 - with gradual adjustment
 - 25723 JETPULL 2 - with gradual adjustment - latex free
 - 25724 JETPULL BABY - with gradual adjustment
- Easy to use, 2 choices for pressure release: instant or gradual release, thanks to the button on the lock. Ribbon in special material to avoid any skin inflammation. Blood proof. Made in Germany.
- 25725 STRIP-VELCRO TOURNIQUET
- Self applicable tourniquet made of a nylon band threads through plastic "D" buckle. Washable and autoclavable.



GIMA code	TOURNIQUETS	Material	Colour	Length cm	Fastening systems	Minimum order	Origin
25720	TUBOLARE Tourniquet	Latex	Yellow	33	Manual	200 pcs.	Europe
25721	GIMA Tourniquet	Rubber-cotton	Yellow	40	Velcro	10 pcs.	Asia
25722	JETPULL 2 Tourniquet	Rubber-cotton	Blue	56	Double lock system: quick and slow release	-	Europe
25723	JETPULL 2 Tourniquet	Rubber-cotton	Green	56	Double lock system: quick and slow release	-	Europe
25724	JETPULL BABY Tourniquet	Rubber-cotton	Child design	56	Double lock system: quick and slow release	-	Europe
25725	STRIP-VELCRO Tourniquet	Velcro	Green	45	Velcro+ring: use with one hand only	-	Asia
25727	FAST Tourniquet	Rubber-cotton	Red	40	Single lock system: quick release	-	Asia
25728	FAST Tourniquet	Rubber-cotton	Green	40	Single lock system: quick release	-	Asia


SYRINGES, HYPODERMIC AND MESOTHERAPY NEEDLES, SCALP VEIN SETS



MADE IN ITALY

STERILE

Made in Italy



STERILE

Made in Italy

Capacity ml	Gauge	3 pieces Syringes with needle*	3 pieces Syringes without needle	Minimum order
0.5	30G 1/3" R0	26000	-	box of 100
1	30G 1/3" R0	26002	-	box of 100
1	25G 5/8"	26005	-	box of 100
2	22G 1 1/4"	-	-	box of 100
2,5	22G 1 1/4"	26010	-	box of 100
5	22G 1 1/4"	26012	-	box of 100
10 CEN	21G 1 1/2"	26015	-	box of 80
10 L/L	-	-	26031	box of 100
10 ECC	21G 1 1/2"	26016	-	box of 80
20 ECC	21G 1 1/2"	26020	26032	box of 80
30 L/L	-	-	26034	box of 80
60 CAT	-	-	26036	box of 55
100 CAT	-	-	26038	box of 50

*R0: zero dead space (no residual)
CAT: catheter (large cone)
1/3": lenght in inches of stainless needle (cannula)*
*L/L: Luer Lock
ECC: Eccentric Luer
CEN: Central Luer*

GIMA code	HYPODERMIC NEEDLES - STERILE			Minimum Order
	Gauge	Gauge mm	Colour code	
26050	18G 1 1/2"	1.20x40	pink	box of 100
26051	19G 1 1/2"	1.10x40	cream	box of 100
26052	20G 1 1/2"	0.90x40	yellow	box of 100
26053	21G 1 1/2"	0.80x40	green	box of 100
26054	22G 1 1/4"	0.70x30	black	box of 100
26055	25G 5/8"	0.50x16	orange	box of 100
26056	27G 1/6" MESO	0.40x4	grey	box of 100
26057	27G 1/4" MESO	0.40x6	grey	box of 100

SCALP VEIN SETS (BUTTERFLY TYPE)- STERILE- TUBE 30 CM LONG				
26060	21G L/L	0.80x19	green	box of 100
26061	23G L/L	0.60x19	dark blue	box of 100

CATETHER NEEDLES - STERILE				
26067	20G L/L	32 mm	pink	box of 100
26068	22G L/L	19 mm	dark blue	box of 100

SYRINGES AND NEEDLES: Gima propose a vast range of sterile, latex free syringes with or without needle made in polypropylene and cannula (needles) made in stainless

steel AISI 304. All products are sterilized by Ethylene Oxide 10%+CO₂, 90% method and made in Italy. Box in 15 languages.
* Needle detachable, excluding 26000 and 26002 (fixed)

BD SHARPS CONTAINERS


SHARPS CONTAINERS BECTON-DICKINSON
Space saving containers, autoclavable at 121°C with safety cover. Conforms to BS7320:1990, ADR UN 3291

- 25822 BD SHARPS CONTAINER 1,4l
- 25823 BD SHARPS CONTAINER 3l
- 25827 WALL SUPPORT for 25823
- 25828 BD SHARPS CONTAINER 7l - 19,5x27x25h cm
- 25829 BD SHARPS CONTAINER 22,7l - 32,5x21,5x33h cm



5 ml SYRINGE - RETAIL BOX

- 25983 SYRINGE 5ml
- With needle 23Gx1 1/4" 0.60x32mm. Sterile, packed in single blister. Box of 10 syringes.
Minimum order: carton of 30 boxes of 10 syringes



SHARPS BINS

- 25818 CONTAINER 1.5l - box of 100
 - 25819 CONTAINER 3l - box of 50
- Disposable sharps containers 1,5 and 3 litres



SKIN MARKER

33181 CLINICON - sterile - box of 10
The Clinicon® Surgical Skin Marker features a versatile tip for precise fine line marking preferred in plastic surgery or when held at an angle, broad line marking for applications in general, neuro and vascular surgery, as well as radiology and radiation therapy. The marker's Gentian violet ink is nontoxic, non tattooing and surgical prep-resistant. Markings can be removed from the skin with isopropyl (rubbing) alcohol. The gamma radiation sterilized package is easy to open, HIBC bar coded and includes a separate six inch, flexible ruler calibrated in both inches and centimeters.

