

GIMA: DISPOSABLE S/S SCALPELS & STAINLESS OR CARBON STEEL BLADES

"GIMA" DISPOSABLE BLADES - sterile - dispenser box of 100
 Standard surgical blades are available in high quality carbon steel and stainless steel with or without safety round tip. Manufactured in accordance with the requirements of British Standards BS 2982:1992, BS 5194: 1992 and ISO 7740:1985. Each blade is hallmarked and individually packed in aluminium foil. Carbon Steel blades are VCI paper protected. Blade size are clearly identified on the reverse side of the foil peel-pack and by a colour coded sticker on unit boxes.



27020 Standard blade n. 10

27030 Safety blade n. 10

"GIMA" DISPOSABLE SCALPELS - sterile - box of 10
 S/S blade fitted into an ABS plastic handle. Individually packed in sterilized envelope. Safety blade guard for premium model.



27040-27048


PREMIUM MODEL



27060-27066

STANDARD MODEL

"GIMA" SAFETY DISPOSABLE SCALPELS - sterile
 Scalpels with stainless steel retractable blade with steadfast disposal lock. Supplied in box of 10 pcs.



27050-27053

SAFETY MODEL

SAFETY HANDLE

SHARP INJURY PREVENTION WITH SAFETY HANDLE AND BLADES

- 26915 SAFETY S/S HANDLE N 3 for blades 10-15 Innovative handle allows hand-free removal of blade in compliance with OSHA and CDC
- 26916 SAFETY S/S HANDLE N 4 for blades 20-25




Safety handle N°3 with safety blade N°10


PARAGON BLADES AND SCALPELS

"PARAGON" DISPOSABLE S/S BLADES - sterile
 High quality blades manufactured with best British stainless steel guarantee excellent sharpening and perfect cut control. Individually packed in metal foil. Gamma rays sterilized.
 "PARAGON" DISPOSABLE SCALPELS WITH S/S BLADE - sterile
 Grip on the plastic handle allows stability for long straight incisions. Individually packed in sterilized envelope, reinforced with blade protection, in cardboard box of 10 units.

Blades box of 100	Scalpels box of 10	STERILE, S/S BLADE
26901	26920	Blade 10
26902	26921	Blade 11
26903	26922	Blade 12
26904	26923	Blade 15
26905	26924	Blade 20
26906	26925	Blade 21
26907	26926	Blade 22
26908	26927	Blade 23
26909	26928	Blade 24



26920-26928



GIMA code	SCALPEL HANDLES - STAINLESS STEEL			Minimum order
26913	Scalpel Handle N° 3 for blades 10-15			10 pcs
26914	Scalpel Handle N° 4 for blades 20-25			10 pcs
Carbon steel box of 100	Stainless steel box of 100	Safety stain. steel box of 100	GIMA STERILIZED BLADES	
27000	27020	27030	Sterilized blades No.10	
27001	27021	-	Sterilized blades No.11	
27002	27022	-	Sterilized blades No.12	
27003	27023	27033	Sterilized blades No.15	
27009	-	-	Sterilized blades No.15c	
27004	27024	27034	Sterilized blades No.20	
27005	27025	27035	Sterilized blades No.21	
27006	27026	27036	Sterilized blades No.22	
27007	27027	27037	Sterilized blades No.23	
27008	27028	27038	Sterilized blades No.24	
Premium Model	Standard Model	Safety Model	GIMA STERILIZED DISPOSABLE SCALPELS STAINLESS STEEL (PACK OF 10 PCS.)	
27040	27060	27050	Disposable scalpels - blade 10	
27041	27061	27051	Disposable scalpels - blade 11	
27042	-	-	Disposable scalpels - blade 12	
27043	27063	27053	Disposable scalpels - blade 15	
27049	-	-	Disposable scalpels - blade 15c	
27044	27064	-	Disposable scalpels - blade 20	
27045	27065	-	Disposable scalpels - blade 21	
27046	27066	-	Disposable scalpels - blade 22	
27047	27067	-	Disposable scalpels - blade 23	
27048	27068	-	Disposable scalpels - blade 24	

SURGICAL PREP USE RAZORS



GIMA code	NON STERILE DISPOSABLE RAZORS stainless steel blades	Minimum order
27080	Single blade - 102x40.3 mm	box of 100
27081	Double blade - 100x40.5 mm	box of 100
27082	Single blade with comb	box of 50

KAI BIOPSI BLADE



• 26935 BIOPSI BLADE - box of 20 pcs
 A new biopsy product from the experience of KAI: blade with sure grip handle for precise biopsy. Sterile and ready to use

SUITABLE FOR SAUCERIZATION

SINGLE-USE STERILE DERMAL CURETTES



Stainless steel semi-sharp cutting edge for abnormal soft or weak tissue. Size embossed on each curette for easy identification. Sterile ready-to-use, use once and discard. Eliminates need for re-sharpening and sterilization. Wide range of sizes. **Used in the following minor surgical skin procedures:** Curettage of a basal cell carcinoma, curettage of a verruca and curettage of any other benign neoplasms of the skin. **Clinical Areas.** Dermatology, surgical, podiatry, veterinary

GIMA code	DERMAL CURETTES	Minimum order
26950	Dermal Curettes Ø 2 mm	20 pcs.
26951	Dermal Curettes Ø 3 mm	20 pcs.
26952	Dermal Curettes Ø 4 mm	20 pcs.
26953	Dermal Curettes Ø 5 mm	20 pcs.
26954	Dermal Curettes Ø 7 mm	20 pcs.




SINGLE-USE STERILE BIOPSY PUNCHES - LATEX FREE

Stainless steel sharp cutting edge, ribbed handle for sure grip. Size embossed on each punch for easy identification. Suggested for dermatology, gynaecology, podiatry. Each punch is individually sealed in a peel-apart see-through package.

*Internal Plunger System
Ejects lodged skin specimen inside
the metal lumen of the punch*



Standard	With plunger	BIOPSY PUNCHES	Blade length mm	Minimum order
26939	26970	Biopsy Punches Ø 1 mm	7	20 pcs.
26940	26971	Biopsy Punches Ø 1.5 mm	7	20 pcs.
26941	26972	Biopsy Punches Ø 2 mm	7	20 pcs.
26938	-	Biopsy Punches Ø 2.5 mm	7	20 pcs.
26942	26974	Biopsy Punches Ø 3 mm	7	20 pcs.
26943	-	Biopsy Punches Ø 3.5 mm	7	20 pcs.
26944	26976	Biopsy Punches Ø 4 mm	7	20 pcs.
26946	-	Biopsy Punches Ø 5 mm	7	20 pcs.
26947	-	Biopsy Punches Ø 6 mm	7	20 pcs.
26948	-	Biopsy Punches Ø 8 mm	7	20 pcs.



*Seamless Stainless Steel
Cutting Edge
Ensures a consistently
smooth cutting surface,
allowing optimal speci-
men quality*




SUTURES - FINE WOUND CLOSURE MATERIAL

CLINISILK STERILE BLACK BRAIDED SILK SUTURE B.P.
Non-absorbable sutures made from pure natural silk, braided to form a multifilament thread, double treated with a silicone coating to prevent capillarity and to enhance ease of passage through tissue. Silk sutures are considered to have the best handling and tying characteristics. Suture length: 45 cm. Application: soft tissue approximation and ligation.

CLININY STERILE BLUE MONOFILAMENT NYLON SUTURE B.P.
Non-absorbable synthetic monofilament suture made from Polyamide 6.6. No tissue reaction, excellent tensile strength, easy handling and secure knotting. Suture length: 45 cm. Nylon is a highly favoured material for skin closure. Application: skin closure, general soft tissue approximation and ligation.

CLINISORB STERILE VIOLET BRAIDED SYNTHETIC ABSORBABLE SUTURE
CliniSorb is a synthetic absorbable braided suture made from a polymer of Glycolic acid. This material is inert, non collagenous, non-antigenic, non-pyrogenic and has a high tensile strength and predictable absorption rate, allowing wound support for approximately 21 days. It is treated with a surface coating to ensure smooth tissue passage and enhance knot sliding. It is not affected by infection, inflammatory processes or exposure to protein digesting enzymes in the gastro intestinal tract. Suitable for both skin closure and where non-permanent wound support is indicated. Application: dental, general skin closure gynaecology, plastic and general surgery, orthopedics and urology.



GIMA code		CLINISILK AND CLININY SUTURES - pack of 12 pcs.		
Silk	Nylon	Gauge	Needle	Type of needle
25910	25920	2/0	26 mm	Reverse cutting - 3/8
25911	-	3/0	22 mm	Conventional cutting - 1/2
-	25921	3/0	20 mm	Conventional cutting - 3/8
25912	25922	3/0	24 mm	Reverse cutting - 3/8
-	25923	3/0	25 mm	Fine conventional cutting - 3/8
25914	25924	4/0	19 mm	Reverse cutting - 3/8
-	25925	5/0	16 mm	Reverse cutting - 3/8
25916	25926	5/0	19 mm	Reverse cutting - 3/8
-	25927	6/0	12 mm	Reverse cutting - 3/8
25918	25928	6/0	16 mm	Reverse cutting - 3/8

GIMA code		CLINISORB SUTURES - pack of 12 pcs. - synthetic		
Gauge	Needle	Type of needle		
25900	2/0	24 mm	Reverse cutting - 3/8 - 70 cm	
25901	3/0	19 mm	Reverse cutting - 3/8 - 45 cm	
25902	4/0	19 mm	Reverse cutting - 3/8 - 45 cm	

NELATON CATHETERS

Male: two lateral eyes with smooth edges

Female: two lateral eyes with smooth edges

GIMA code	DEHP FREE STERILE NELATON CATHETER two lateral eyes with smooth edges Single use	Minimum order
MALE - length 38 cm		
22003	Male nelaton catheter ch/fr 10	box of 100
22004	Male nelaton catheter ch/fr 12	box of 100
22005	Male nelaton catheter ch/fr 14	box of 100
22006	Male nelaton catheter ch/fr 16	box of 100
FEMALE - length 18 cm		
22023	Female nelaton catheter ch/fr 10	box of 100
22024	Female nelaton catheter ch/fr 12	box of 100
22025	Female nelaton catheter ch/fr 14	box of 100
22026	Female nelaton catheter ch/fr 16	box of 100

FOLEY CATHETERS

Three eyes with smooth edges

Two simmetrical eyes with smooth edges

GIMA code	STERILE 100% SILICONE COATED FOLEY CATHETERS - SINGLE USE	Minimum order
Two-way catheters - baloon 5 or 30 ml - soft valve		
22063	2-way foley catheter ch/fr 10 - baloon 5 ml	box of 50
22064	2-way foley catheter ch/fr 12 - baloon 30 ml	box of 50
22065	2-way foley catheter ch/fr 14 - baloon 30 ml	box of 50
22066	2-way foley catheter ch/fr 16 - baloon 30 ml	box of 50
22067	2-way foley catheter ch/fr 18 - baloon 30 ml	box of 50
Three-way catheters - baloon 30 ml - soft valve		
22078	3-way foley catheter ch/fr 20	box of 50
22079	3-way foley catheter ch/fr 22	box of 50
22080	3-way foley catheter ch/fr 24	box of 50

SUCTION CATHETERS

Length 90 cm

Length 50 cm	DEHP FREE STERILE SUCTION CATHETER with conical adapter	Minimum order
22110	22120 Suction catheter ch/fr 6	box of 100
22112	22122 Suction catheter ch/fr 8	box of 100
22113	22123 Suction catheter ch/fr 10	box of 100
22114	22124 Suction catheter ch/fr 12	box of 100
22115	22125 Suction catheter ch/fr 14	box of 100

FEEDING TUBES

Marked between 15-20 cm from the distal end

GIMA code	DEHP FREE STERILE FEEDING TUBE length: 50 cm - with luer fitting and sealing plug adapter	Minimum order
22150	Feeding tube ch/fr 5	box of 100
22151	Feeding tube ch/fr 6	box of 100
22152	Feeding tube ch/fr 8	box of 100
22153	Feeding tube ch/fr 10	box of 100

URINE BAGS MADE IN EUROPE

Internal Ø 4.8 mm

GIMA code	URINE BAGS	Minimum order
28700	28700 Urine bag 2000 cc, tube 90 cm	box of 180
28701	28701 Urine bag 2000 cc, tube 90 cm with non-return/exhaust valve 18 pouches of 10 pcs.	box of 180
28702	28702 Urine bag 2000 cc, tube 130 cm	box of 180
28703	28703 Urine bag 2000 cc, tube 130 cm with non-return/exhaust valve 18 pouches of 10 pcs.	box of 180
28706	28706 Urine leg bag 750 cc, tube 35 cm complete with lacets and buttons - 30 pouches of 10 pcs.	box of 300
28710	28710 Men's external catheter Ø 3 mm - latex	box of 30
28715	28715 Biadhesive tape 2x15x0,5 cm	box of 30

BIOLOGICAL COLLECTORS

PERSONALIZED PACKAGING FOR 20,000 PCS.

GIMA code	POLYPROPYLENE CONTAINERS WITH SINGLE BAG, SINGLE BOX AND ADHESIVE LABEL	Minimum order
25970	25970 Urine container - 120 ml with screw cap - graduated to 100 ml	box of 120
25971	25971 Biological products container - 30 ml screw cap with small spatula	box of 100
GIMA code	POLYSTYRENE CONTAINERS WITH SINGLE BAG, SINGLE BOX AND ADHESIVE LABEL	Minimum order
25972	25972 Urine test tube - 12 ml Blue sanitary strip	box of 120

TRAVEL WC

- 25985 PRONTO TRAVEL WC - box of 3 pcs
- Practical bag that instantly gel patient urine and absorb gel
- No closure needed
- Waste disposal safe, simply throw away
- Reusable until it is full
- Super absorbent polymer pouch included



TRAVELLING ENEMA

28720	• 28720 TRAVELLING ENEMA Supplied in carton box with cannula and tubes. Made in Italy.
-------	--

24 HOURS URINE CONTAINER

- 25974 24 HOURS URINE CONTAINER
Graduated to 2,500 ml - box of 24
Male-female container.
Each container is single packed in carton box.




ROUND STERILIZER DRUMS - STAINLESS STEEL AISI304

ROUND STERILIZER DRUM for autoclave and hot air sterilizer. Sliding removable band and drop down clasps.

GIMA code	STERILIZER DRUMS
26673	Ø 120x h 120 mm
26674	Ø 140x h 140 mm
26676	Ø 160x h 120 mm
26678	Ø 180x h 180 mm
26680	Ø 190x h 160 mm
26681	Ø 240x h 160 mm
26682	Ø 240x h 240 mm
26683	Ø 290x h 160 mm
26684	Ø 290x h 240 mm
26687	Ø 290x h 290 mm
26688	Ø 340x h 240 mm
26685	Ø 340x h 340 mm
26686	Ø 390x h 250 mm



STERILIZATION CONTAINER SYSTEM



STANDARD CONTAINER WITH FILTER*			
LID COLOUR	Grey	Red	Blue
Height mm	[NON PERFORATED BOTTOM]		[PERFORATED BOTTOM*]
		SMALL 285x280 mm	
	1 filter		2 filter
100	37000	37010	37030
135	37001	-	37031
150	37002	37012	37032
200	37003	-	37033
260	37004	-	-
	MEDIUM 465x280 mm		
	1 filter		2 filter
100	37050	37060	37080
135	37051	-	37081
150	37052	37062	37082
	LARGE 580x280 mm		
	2 filter		4 filter
100	37100	37110	37130
135	37101	-	37131
150	37102	37112	37132
200	37103	-	37133
260	37104	-	-

*Perforated bottom allows much quicker sterilization process (half time than with not perforated)

BARRIER SYSTEM CONTAINER WITH VALVE**			
LID COLOUR	Grey	Red	Blue
Height mm	[NON PERFORATED BOTTOM]		[PERFORATED BOTTOM*]
		SMALL 285x280 mm	
100	37150	37160	37180
135	37151	-	37181
150	37152	37162	37182
200	37153	-	37183
260	37154	-	-
	MEDIUM 465x280 mm		
100	37200	37210	37240
135	37201	-	37241
150	37202	37212	37242
	LARGE 580x280 mm		
100	37250	37260	37280
135	37251	-	37281
150	37252	37262	37282
200	37253	-	37283
260	37254	-	-

**This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.




WIRE BASKET			
Size mm	Height		
	50mm*	70mm**	100mm***
255x245	37301	37302	37303
405x255	37306	37307	37308
540x255	37311	37312	37313

*Suggested for container with height 100/135 mm
 ** Suggested for container with height 135/150 mm
 *** Suggested for container with height >150 mm



GIMA code	ACCESSORIES
37320	Silicone mat 220x230 mm - perforated
37321	Silicone mat 380x230 mm - perforated
37323	Silicone mat 520x230 mm - perforated
37325	Filter holder for paper
37326	Filter holder for textile
37328	Paper filter - box of 100 - for standard container
37329	Textile filter - for standard container
37330	Plastic security seal - box of 100
37331	Paper label with indicator - box of 100

ALUMINIUM BOXES



- 26662 ALUMINIUM BOX 17.5 x 7.6 x h 2 cm
 - 26663 ALUMINIUM BOX 18.5 x 9.5 x h 3 cm
 - 26664 ALUMINIUM BOX 21.8 x 10.6 x h 3 cm
 - 26665 ALUMINIUM BOX 21.8 x 10.6 x h 5 cm
- Useful containers for the sterilization of surgical instruments and to maintain them sterile. Rubber tray not included. Nice aluminium box with blue base and silver cover.
- 26672 RUBBER TRAY 20x10 cm (for 26669, 26664, 26665)

STAINLESS STEEL BOXES



- STAINLESS STEEL BOXES**
- 26654 STAINLESS STEEL BOX 20 x 10 x h 6 cm
 - 26655 STAINLESS STEEL BOX 25 x 12 x h 6 cm
 - 26656 STAINLESS STEEL BOX 30 x 15 x h 6 cm
- STAINLESS STEEL BOXES WITH HANDLES**
- 26668 STAINLESS STEEL BOX 18 x 8 x h 4 cm
 - 26669 STAINLESS STEEL BOX 20 x 10 x h 4.5 cm
 - 26670 STAINLESS STEEL BOX 25 x 12.5 x h 4.6 cm
 - 26671 STAINLESS STEEL BOX 30 x 15 x h 6 cm
- Autoclavable 100°

PLASTIC HOLLOWARE



AUTOCLAVABLE AT 135°C


PLASTIC HOLLOWARE - autoclavable hundreds of time
For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% cheaper than stainless steel equivalent
- autoclavable at 135°C with no limits
- attractive modern design
- lighter than equivalent stainless steel

All products are manufactured from polypropylene material That is easy to clean, resistant to most healthcare chemicals.

- 26630 GALLIPOT/LOTION BOWL 60 mm - 98 ml
- 26631 GALLIPOT/LOTION BOWL 80 mm - 200 ml
- 26632 GALLIPOT/LOTION BOWL 100 mm - 350 ml
- 26633 GALLIPOT/LOTION BOWL 150 mm - 1000 ml
- 26635 KIDNEY DISH 200 x 98 mm - 500 ml
- 26636 KIDNEY DISH 250 x 130 mm - 750 ml with special nodules on rim to prevent vacuum during autoclaving
- 26638 INSTRUMENT TRAY 200 x 150 x 51 mm
- 26639 INSTRUMENT TRAY 300 x 250 x 52 mm
- 26640 COMPARTMENT TRAY 270 x 180 x 41 mm



26640

DISPOSABLE KIDNEY DISH



- 26645 DISPOSABLE KIDNEY DISH - box of 300 pcs.
- Safe, clean, easy to use and disposal. This kidney dish can be used for any kind of medical purpose.

26645

TECHNICAL SPECIFICATIONS

Material:	molded paper pulp, waterproof
Colour:	typical grey-brown
Packaging:	masters cartons of 300 pieces (3 bags of 100 pcs.)
Size:	250x136x45 mm
Maximum capacity:	about 700 ml
Disposal:	biodegradable after being discarded, in line with environmental protection requirements marked CE - Class I non sterile
Medical device:	

PLASTIC BED PAN - URINAL



25979



25978

- 25979 PLASTIC BED PAN - box of 20
- 25978 PLASTIC MEN URINAL - box of 60

SURGERY COTTON DRAPES



- 26495 DRAPE 90X150
 - 26496 DRAPE 150x150
 - 26497 DRAPE 250x150
- Made in high quality 100% cotton, green colour. Made in Italy.

STAINLESS STEEL HOLLOWARE - AUTOCLAVABLE 121°C - AISI 304


GIMA code	INSTRUMENT TRAYS
Without cover	
26600	208x109x h 15 mm
26601	210x160x h 25 mm
26599	223x126x h 45 mm
26602	264x172x h 47 mm
26603	306x196x h 50 mm
26604	355x254x h 50 mm
26597	380x304x h 50 mm
26598	440x320x h 64 mm

GIMA code	LOTION BOWLS
Round Bowls with lip	
26616	Ø 56xh 32 mm - 50 ml
26617	Ø 88xh 46 mm - 170 ml
26618	Ø 128xh 60 mm - 450 ml
26619	Ø 158xh 75 mm - 940 ml
26621	Ø 208xh 84 mm - 1940 ml
26622	Ø 258xh 92 mm - 3200 ml



GIMA code	With cover
26605	223x126x h 45 mm
26606	264x172x h 47 mm
26607	306x196x h 50 mm
26608	355x254x h 50 mm
26609	440x320x h 64 mm

GIMA code	DRESSING JAR WITH LID
26623	0.5 l Ø 102x65 mm
26624	1 l Ø 103x126 mm
26625	2 l Ø 127x161 mm



GIMA code	BED PANS AND URINALS
Bed Pan	
26592	320x85 mm (Round with cover)
26593	354x289x100 mm (Perfection Type)
26590	418x292x85 mm with ring (Fracture Type)
26594	418x292x85 mm with handle (Fracture type)
Urinals	
26591	Ø 114x220 mm - 1.5l (Male Vertical)
26578	Male urinal - 1.2l
26579	Female urinal - 1.2l

GIMA code	THERMOMETER JAR
26628	Ø 33x h 80 mm

GIMA code	FORCEPS JAR
26626	Ø 55x h 140 mm
26627	Ø 55x h 180 mm



GIMA code	KIDNEY DISHES
Shallow	
26610	207x128x h 33 mm
26611	254x141x h 33 mm
26612	280x141x h 33 mm
Deep	
26580	162x77x31 mm
26581	207x98x39 mm
26582	247x122x43 mm
26583	309x149x59 mm
Deep with cover	
26584	207x98x39 mm
26585	247x122x43 mm
26589	309x149x59 mm



































GIMA code	WASH BASINS WITHOUT LIP
26586	Ø 310x h 72 mm - 3l
26587	Ø 318x h 84 mm - 4l
26588	Ø 405x h 95 mm - 7l

TUNGSTEN CARBIDE INSTRUMENTS - GOLD LINE







GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments to help surgeons to perform the precision surgeries with best results. These instruments are more reliable for hardness and corrosion resistance.











NEEDLE HOLDERS	SCISSORS	WIRE CUTTERS
 26544 Webster needle holder - 12.5 cm	 26500 Straight scissors blunt/sharp - 14.5 cm 26501 Straight scissors blunt/blunt - 14.5 cm 26502 Straight scissors sharp/sharp - 14.5 cm  26503 Curved scissors blunt/sharp - 14.5 cm 26504 Curved scissors blunt/blunt - 14.5 cm 26505 Curved scissors sharp/sharp - 14.5 cm	 26570 Wire Cutter - 14 cm - for soft wires 0-1 mm  26571 Wire Cutter - 18 cm - for hard wires cutting up to 1.6 mm 26572 Wire Cutter - 23 cm - for hard wires cutting up to 2.0 mm  26575 Wire holding - 18 cm
 26546 Crile Wood needle holder - 15 cm	 26526 Iris scissors - straight - 11.5 cm  26527 Iris scissors - curved - 11.5 cm	<h3>ORGAN HOLDING FORCEPS</h3>  26564 Duval - 20 cm
 26530 Mayo Hegar needle holder - 16 cm  26531 Mayo Hegar needle holder - 18 cm  26532 Mayo Hegar needle holder - 20 cm	 26528 Goldman Fox scissors - straight - 13 cm  26529 Goldman Fox scissors - curved - 13 cm	<h3>MICRO NEEDLE HOLDER</h3>  Plane tip Rough tip 26550 Castroviejo Needle Holder straight - 14 cm - plane tip 26552 Castroviejo Needle Holder curved - 14 cm - plane tip 26551 Castroviejo Needle Holder straight - 18 cm - rough tip
 26540 Olsen Hegar needle holder - 14 cm  26541 Olsen Hegar needle holder - 16.5 cm	 26510 Metzenbaum scissors - straight - 14 cm  26511 Metzenbaum scissors - straight - 18 cm  26512 Metzenbaum scissors - straight - 20 cm  26514 Metzenbaum scissors - curved - 14 cm  26515 Metzenbaum scissors - curved - 18 cm  26516 Metzenbaum scissors - curved - 20 cm	<h3>ADSON FORCEPS</h3>  26554 Castroviejo - 13 cm 26560 Adson - 12 cm 26561 Adson 1x2 teeth - 12 cm
 26536 Mathieu needle holder - 14 cm  26537 Mathieu needle holder - 17 cm	 26520 Mayo Scissors - straight - 14.5 cm  26521 Mayo Scissors - straight - 18 cm  26522 Mayo Scissors - curved - 14.5 cm  26523 Mayo Scissors - curved - 18 cm	

SCALPELS

SCALPELS AND PROBES
 26698 Scalpel straight blade - 13 cm  26703 Scalpel straight blade - 17.5 cm
 26699 Scalpel pot-bellied blade - 13 cm  26704 Scalpel pot-bellied blade - 17.5 cm
 26701 Probe - 15 cm - Ø 2 mm
 26702 Butterfly probe - 14 cm

NEEDLE HOLDERS

NEEDLE HOLDERS
 26852 Crile Wood needle holder - 15 cm
 26818 Mayo Hegar needle holder - 14 cm  26732 Mayo Hegar needle holder - 16 cm  26819 Mayo Hegar needle holder - 18 cm  26733 Mayo Hegar needle holder - 20 cm  26741 Mathieu needle holder - 14 cm  26730 Mathieu needle holder - 16 cm  26731 Mathieu needle holder - 20 cm

MEDICATION

STITCH REMOVER FORCEPS/ SCISSORS
 26737 Spencer scissors - 13 cm for suture
 26722 Michel stitch remover forceps - 13 cm
LISTER BANDAGE SCISSORS
 34125 Lister scissors - 11 cm 34129 Lister scissors - 15 cm 34130 Lister scissors - 18 cm

STAINLESS STEEL SURGICAL INSTRUMENTS

TOWEL FORCEPS



- 26743 Backhaus forceps - 9 cm
- 26736 Backhaus forceps - 11 cm
- 26747 Backhaus forceps - 13 cm



- 26739 Backhaus Forceps for TNT- 13 cm

DISSECTING FORCEPS



- 26690 Anatomy forceps - 12 cm
- 26705 Anatomy forceps - 14 cm
- 26691 Anatomy forceps - 16 cm
- 26706 Anatomy forceps - 18 cm
- 26752 Anatomy forceps - 20 cm




- 26692 Surgery forceps - 12 cm - 1x2
- 26707 Surgery forceps - 14 cm - 1x2
- 26693 Surgery forceps - 16 cm - 1x2
- 26708 Surgery forceps - 18 cm - 1x2
- 26753 Surgery forceps - 20 cm - 1x2




- 26750 Adson forceps - 12 cm
- 26751 Adson forceps - 12 cm - 1x2
- 26694 Micro Adson forceps - 12 cm
- 26695 Micro Adson forceps - 12 cm - 1x2

DRESSING SPONGE HOLDING FORCEPS



- 26745 Pean forceps - 14 cm
- 26720 Pean forceps - 16 cm
- 26710 Pean forceps - 20 cm



- 26800 Foerster forceps - 20 cm
- 26801 Foerster forceps - 25 cm



- 26805 Foerster "small ring" forceps - 25 cm



- 26845 Duval 14 cm x 10 mm

SCISSORS



Curved

Straight Blunt/Sharp

- 26754 Straight scissors blunt/sharp - 11.5 cm
- 26725 Straight scissors blunt/sharp - 14.5 cm
- 26729 Straight scissors blunt/sharp - 16 cm
- 26853 Straight scissors blunt/sharp - 18 cm
- 26726 Straight scissors blunt/sharp - 20 cm
- 26755 Curved scissors blunt/sharp - 11.5 cm
- 26727 Curved scissors blunt/sharp - 14.5 cm
- 26756 Curved scissors blunt/sharp - 16 cm
- 26854 Curved scissors blunt/sharp - 18 cm
- 26728 Curved scissors blunt/sharp - 20 cm




Straight Blunt/Blunt

- 26855 Straight scissors - blunt/blunt - 11.5 cm
- 26740 Straight scissors - blunt/blunt - 14.5 cm
- 26856 Straight scissors - blunt/blunt - 16 cm
- 26857 Straight scissors - blunt/blunt - 18 cm
- 26887 Straight scissors - blunt/blunt - 20 cm
- 26888 Curved scissors - blunt/blunt - 11.5 cm
- 26742 Curved scissors - blunt/blunt - 14.5 cm
- 26889 Curved scissors - blunt/blunt - 16 cm
- 26890 Curved scissors - blunt/blunt - 18 cm
- 26891 Curved scissors - blunt/blunt - 20 cm



Straight Sharp/Sharp

- 26892 Straight scissors - sharp/sharp - 11.5 cm
- 26744 Straight scissors - sharp/sharp - 14.5 cm
- 26893 Straight scissors - sharp/sharp - 16 cm
- 26894 Straight scissor - sharp/sharp - 18 cm
- 26895 Straight scissor - sharp/sharp - 20 cm
- 26896 Curved scissor - sharp/sharp - 11.5 cm
- 26746 Curved scissor - sharp/sharp - 14.5 cm
- 26897 Curved scissor - sharp/sharp - 16 cm
- 26898 Curved scissor - sharp/sharp - 18 cm
- 26899 Curved scissor - sharp/sharp - 20 cm



- 26822 Iris scissors - straight - 11 cm
- 26824 Iris scissors - curved - 11 cm

SURGICAL SCISSORS




- 26840 Metzenbaum scissors - straight 14.5 cm
- 26841 Metzenbaum scissors - straight 18 cm
- 26843 Metzenbaum scissors - curved 14.5 cm
- 26844 Metzenbaum scissors - curved 18 cm



- 26846 Mayo-stille scissors - straight 14.5 cm
- 26847 Mayo-stille scissors - straight 18 cm
- 26848 Mayo-stille scissors - straight 20 cm
- 26849 Mayo-stille scissors - curved 14.5 cm
- 26850 Mayo-stille scissors - curved 18 cm
- 26851 Mayo-stille scissors - curved 20 cm

HAEMOSTATIC FORCEPS



- 26714 Mosquito forceps - straight - 12.5 cm
- 26872 Mosquito forceps - straight - 14 cm
- 26715 Mosquito forceps - straight - 16 cm
- 26716 Mosquito forceps - curved - 12.5 cm
- 26873 Mosquito forceps - curved - 14 cm



- 26757 Kelly forceps - straight - 14 cm.
- 26717 Kelly forceps - straight - 16 cm.
- 26758 Kelly forceps - curved - 14 cm.
- 26718 Kelly forceps - curved - 16 cm.



- 26828 Kocher forceps - straight - 14 cm - 1x2
- 26711 Kocher forceps - straight - 16 cm - 1x2
- 26829 Kocher forceps - straight - 18 cm - 1x2
- 26712 Kocher forceps - straight - 20 cm - 1x2
- 26831 Kocher forceps - curved - 14 cm - 1x2
- 26832 Kocher forceps - curved - 16 cm - 1x2
- 26833 Kocher forceps - curved - 18 cm - 1x2
- 26834 Kocher forceps - curved - 20 cm - 1x2



- 26749 Klemmer forceps - straight - 14 cm
- 26709 Klemmer forceps - straight - 16 cm
- 26835 Klemmer forceps - straight - 18 cm
- 26713 Klemmer forceps - straight - 20 cm
- 26836 Klemmer forceps - curved - 14 cm
- 26837 Klemmer forceps - curved - 16 cm
- 26838 Klemmer forceps - curved - 18 cm
- 26839 Klemmer forceps - curved - 20 cm



- 26734 Crile forceps - straight - 16 cm
- 26735 Crile forceps - curved - 16 cm




- 26721 Allis forceps - 15 cm
- 26759 Allis forceps - 19 cm

MAGILL AND COLLIN FORCEPS
see at page 222



GIMA code	INSTRUMENTS KIT IN ALUMINIUM BOX OR CORDURA BAG													
		26701	26702	26703	26704	26705	26707	26709	26722	26723	26724	26725	26730	26732
26760	Trousse standard - cordura bag 9 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26761	Trousse classica - cordura bag 10 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26762	Trousse suprema - cordura bag 11 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26768	Trousse suprema - aluminium box 11 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•


STAINLESS STEEL RETRACTORS



- 26858** Weitlaner-Loktite retractor - sharp
13 cm - 3x4 prongs
- 26859** Weitlaner-Loktite retractor - blunt
13 cm - 3x4 prongs
- 26865** Weitlaner Retractor - sharp
20 cm - 3x4 prongs
- 26866** Weitlaner Retractor - blunt
20 cm - 3x4 prongs



- 26867** Farabeuf retractors - set of 2 - 12 cm




- 26868** Frazier Hook - sharp - 13 cm
- 26869** Frazier Hook - blunt - 13 cm



- 26723** Hook - 15 cm.



- 26860** Masing Retractor - 14 cm



- 26877** Roux Retractor set of 5 (3+2)



- 26862** Volkman- 3 teeth - sharp - 22 cm
- 26863** Volkman- 3 t eeth - blunt - 22 cm



- 26861** Langenbeck Retractor - 21 cm

DERMATOLOGY INSTRUMENTS AND KIT



- 26707** ● Forceps surgery - 14 cm



- 26724** ● Forceps thin end - 12 cm



- 26821** ● Unna round - 14 cm



- 26820** ● Unna Oval - 14 cm



- 26823** ● Round forceps - 9 cm



- 26822** ● Scissors IRIS straight- 11 cm
- 26824** Scissors IRIS curved - 11 cm



- 26825** Katsch - 14 cm




- 26826** Gillies hook - 15 cm - small

- **26830** DERMATOLOGY KIT
- Includes 6 instruments marked with ●




26830

GYNAECOLOGY INSTRUMENTS AND KIT



- 26807** Novak suction cannula - 23 cm




- 26804** ★ Sims histerometer - 33 cm - silver



- 26803** ★ Bosch/Sims scissors - 23 cm




- 26802** ★ Schroeder forceps - 25 cm
(Pozzi-Palmer)



- 26801** ★ Foerster forceps - 25 cm

- **26810** GYNAECOLOGY KIT
- Includes 4 instruments marked with ★ in a nice nylon bag.



- 26813** Sims uterine curette - sharp - 8 mm - 26 cm
- 26814** Sims uterine curette - sharp - 9 mm - 26 cm
- 26815** Sims uterine curette - sharp - 14 mm - 26 cm



- 26817** Set of chrome plated Hegar Dilators
- Include 14 dilators Ø 4-17 mm



26810



- 26806** Schumacher biopsy forceps - 24 cm



- 26808** Schubert biopsy forceps - 21 cm



- 26809** Kogan endospeculum - 24 cm



- 26811** Bozeman dressing forceps straight - 25 cm
- 26812** Bozeman dressing forceps curved - 25 cm

BEAUTY STAINLESS STEEL INSTRUMENTS



- 26870** Cuticle Scissor - 9 cm

- 26871** Nail Scissor - 9 cm

- 26875** Tweezer - 8 cm




- 26885** Cuticle Nipper - 11.5 cm



- 26883** Ingrowing Nail Cutter - 11.5 cm



- 26880** Nail Cutter - 14 cm



- 26881** Nail Cutter Slant Jaws - 14 cm



- 26886** Nail Nipper - 12.5 cm



- 26882** Nail Cutter with 4 hinges - 16 cm

- 26876** Tweezer - 9.5 cm

DENTAL STAINLESS STEEL INSTRUMENTS

IMPRESSION TRAY





60040
set of 10 perforated
impression trays
(5 lower + 5 upper)




60045
set of 10 solid
impression trays
(5 lower + 5 upper)

EXTRACTING-FORCEPS - FOR JAWS



EXPLORER AND PROBES




60101 Explorer - fig. 17/23

60100 Explorer probe

60105 WHO probe

GRACEY AND Mc CALL CURETTES



60120 Mc Call curettes - fig. 13/14


60125 Gracey curette - fig. 1/2 ant.

60127 Gracey curette - fig. 7/8 medium

60130 Gracey curette - fig. 11/12 medium

60131 Gracey curette - fig. 13/14 post.

SCALER



60140 Scaler - fig. 1

60141 Mc Call scaler - fig. 18

60142 Mitchel scaler

60144 Langer scaler - fig. 1/2

60145 Langer scaler - fig. 3/4


60148 Hygienist scaler - fig. H6/7

60150 Scaler - fig. 2045

60155 Periodontia instr. - fig. ST 2/3

60156 Periodontia instr.

CARVER



60170 Cement spatula

60171 Wiland modeling carver

60172 Interproximal modelling carver

60173 Beale modelling carver - mod. 7A

60174 Fahrenstock - modelling carver - 13 cm

FORCEPS




60190 London college forceps - 15 cm

CARTRIDGE SYRINGE



60195 Cartridge Syringe - 2.2 cc
with thread in millimeter
with aspiration



60201

60200 Bein
root
elevator

60202

Seldin
root
elevator

60203

60204 **60205**

Winter
root
elevator

60210

Amalgam
carrier

Molt
Mouth gag
60220 - 11 cm
60221 - 14 cm

60225

Bone Cutter
- 15 cm

60226

Beyer Forceps
- 18 cm

ENT STAINLESS STEEL INSTRUMENTS AND KIT



Tongue depressors

- 25620 Mayo tongue depressor - 17 cm
- 25621 Bruenings tongue depr - fenestrated - 19 cm ●
- 25623 Tobold tongue depressor




Nose/ear speculum

- 26774 Voltolini nose speculum - adjustable - Size 2 ●
- 26775 Voltolini nose speculum - adjustable - Size 0
- 26770 Set of 3 ear speculum Hartmann - Ø 4, 5, 6 mm ●
- 26784 Thudicum nose clip



Nose Ear forceps and scissors

- 26773 Lucae ear forceps - 14 cm ●
- 26783 Troeltsch ear forceps - 12 cm
- 26779 Hartmann ear forceps - 14 cm
- 26781 Heymann scissors - 18 cm



Suction cannulas

- 26795 Frazier nose/ear suction cannula Ø 2 mm (18 cm)
- 26796 Frazier nose/ear suction cannula Ø 3 mm (18 cm)
- 26793 Frazier nose suction cannula Ø 4 mm (18 cm)
- 26794 Frazier nose suction cannula Ø 5 mm (18 cm)
- 26797 Rosen ear suction cannula Ø 1.5 mm (55 mm)
- 26798 Rosen ear suction cannula Ø 3 mm (55 mm)
- 26799 Suction adaptor with LUER tip for 26797/8



Hooks-Loops

- 26771 Lucae Hook - 14 cm (5/4") ●
- 26769 Lucae Hook - fine - 14 cm (5/4")
- 26790 Billeau Loop - 16 cm
- 26785 Swab holder - 16 cm



Killian/Hartmann nose forceps

- 26792 Killian/Hartmann nose forceps - 14 cm (5/4")/15 mm
- 26778 Killian/Hartmann nose forceps - 14 cm (5/4")/22 mm
- 26772 Killian/Hartmann nose forceps - 14 cm (5/4")/35mm ●
- 26776 Killian/Hartmann nose forceps - 14 cm (5/4")/50mm
- 26764 Killian/Hartmann nose forceps - 14 cm (5/4")/75mm
- 26777 F. O. Killian nose forceps - 14 cm/75mm fiber optic
- 26765 Tieck Halle nose speculum - 14 cm


GIMA code mm	Hartmann ear forceps			Picture
	A	B	C	
26782	90	11	0.6	

GIMA code mm	Hartmann micro forceps			Picture
	A	B	C	
26786	80	6	0.3	
26787	80	4	0.6	

● 26788 MICRO EAR CUP-SHAPED POLYPUS FORCEPS
A = 80 mm, B = 4 mm



● 26789 MICRO EAR SCISSORS BELLUCCI
A = 80 mm, B = 4 mm



● 26791 MICRO EAR PUNCH - introduction instrument
A = 80 mm, B = 4 mm



● 26780 ENT KIT Includes all instruments marked with ● inside an aluminium box (Code 26665) 21 x 10 x h 5 cm.

LARYNGEAL MIRRORS AND SPIRIT LAMP



31614



31612

- 31616 ADJUSTABLE MIRROR HANDLE
- 31615 FIXED MIRROR HANDLE

31616

31615



GIMA code	LARYNGEAL MIRROR - STRAIGHT			Minimum order
31600	Mirror number	000	diameter 10 mm	10 pcs. mixed sizes
31601	Mirror number	00	diameter 12 mm	
31602	Mirror number	0	diameter 14 mm	
31603	Mirror number	1	diameter 16 mm	
31604	Mirror number	2	diameter 18 mm	
31605	Mirror number	3	diameter 20 mm	
31606	Mirror number	4	diameter 22 mm	
31607	Mirror number	5	diameter 24 mm	
31608	Mirror number	6	diameter 26 mm	

- 31612 SET 9 LARYNGEAL MIRRORS 000-6
Set of 9 laryngeal mirrors from 000 to 6, all mirrors have a different diameter: from 10 mm to 26 mm.
minimum order: 1 set of 9 mirrors

- 31614 SPIRIT LAMP - 120 cc
Lamp is made of AISI 304 stainless steel.
Brass cap to close flame area to avoid discoloration.

OTOSCILLO ELECTRICAL IRRIGATION DEVICE

• **25815 OTOSCILLO**

Otosillo is designed for automated ear rinsing both in hospitals and during private consultations. Innovative, compact, efficient and rapid, the Otosillo combines the most sophisticated operating qualities with the best security guarantees, all in a single device.

Anywhere, anytime

Otosillo is an electrical irrigation for simplified rinsing of the ear canal. The pulsating stream system makes it easier to remove cerumen, even at low pressure control, Otosillo is particularly adapted to the treatment of sensitive ear canals. It looks less alarming and it is easier to use, with greater accuracy and security, than the alternatives. The water basin is removable to allow easy filling and emptying. The device is readily usable and easy to clean.

Can be mounted on a wall or a cupboard. Made in Germany
Operating power: 220V - 50/60 Hz

Manual in English, French, German, Italian, Spanish

- **25816 SPARE SET (2 tips + disk)**
- **25817 SPARE HANDLE + TUBE**
- **25814 DISPOSABLE TIPS SET**
(100 disposable tips, 1 disk, connector)



25815

GLASS-METAL, ALL-METAL, AND SINGLE-USE EAR SYRINGES



ALL-METAL SCHIMMELBUSCH SYRINGES

High quality, reliable and accurate syringes. Schimmelbusch come complete with 3 cannulas and ear bowl. Individually packed in carton box. Made in Europe.


GLASS-METAL JANET RECORD SYRINGES

High quality, reliable and accurate glass-metal syringes. Janet Record come complete with 2 cannulas, one conical, one oval. Individually packed in carton box. Made in Europe.

GIMA code	ALL-METAL SYRINGES COMPLETE WITH 3 CANNULAS AND EAR BOWL
25806	Schimmelbusch 50 cc
25807	Schimmelbusch 100 cc
25808	Schimmelbusch 200 cc
25809	Spare cannula for Schimmelbush (conical shape)
25812/3	Spare cannula for Schimmelbush (olive shape) small/large

GIMA code	GLASS-METAL SYRINGES COMPLETE WITH 2 CANNULAS
25801	Janet Record 50 cc - 13 x Ø 3.5 cm
25802	Janet Record 100 cc - 14 x Ø 4.7 cm
25803	Janet Record 200 cc
25804	Spare cannula for Janet Record (conical shape)
25805	Spare cannula for Janet Record (olive shape)

HIGH CAPACITY POLYPROPYLENE SYRINGES



26038

Luer adaptor included

HIGH CAPACITY POLYPROPYLENE SYRINGES - sterile

Sterile, single use, latex free syringes for ear and other purposes. Made in Italy.

GIMA code	HIGH CAPACITY, POLYPROPYLENE SYRINGES STERILE	Minimum order
26036	Syringe 60 ml	55 pcs.
26038	Syringe 100 ml	50 pcs.